

# Compression Therapy

Lymphedema • Chronic Edema • Venous Edema



## Our Focus: the Person

Lohmann & Rauscher is focused on the health and care of the individual. Our objective is to provide quality products that aid healthcare providers to make the healing process as comfortable and as short as possible for the patient. This leads to an essential focus on efficient service and product quality.

## Confidence & Partnership

In selecting products from Lohmann & Rauscher, you can rest assured that you have the support and backing of one of the world's leading providers of a variety of medical products, including compression bandaging. We are a reliable partner to meet your customer needs. The breadth and width of our bandaging product line offers many specialized options for a variety of treatment objectives.

## Quality Guarantee

Lohmann & Rauscher's products are assured to be the highest quality available. We guarantee your satisfaction for the normal life of our products.

Please contact our customer service team at **800-279-3863** with any concerns or suggestions.

[www.Lohmann-Rauscher.us](http://www.Lohmann-Rauscher.us)

**Lohmann & Rauscher** is a leading global company that provides quality products across the world within the healthcare industry. Lohmann & Rauscher has a 160-year history as a manufacturer and distributor, specializing in medical devices, compression therapy, advanced wound care and surgical dressings.

**Lohmann & Rauscher International Headquarters are in Rengsdorf, Germany and Vienna, Austria**



Lohmann & Rauscher is the parent company of **L&R USA, Inc.**, with US Headquarters in Milwaukee, Wisconsin, USA. We offer the ease and convenience of "local service" in combination with support from a 160-year global strength. We work as one team to continually focus on innovation, education and research, to meet your needs.

**L&R USA, Inc.**  
**U.S. Headquarters is in Milwaukee, WI**



# Product Listing

## Product (alphabetical listing) Page #

---

|  |    |
|--|----|
| Bandage Roller                                 | 21 |
| Cellona® Synthetic Padding                     | 12 |
| Compression Bandaging for Lymphedema Treatment | 24 |
| Contributions towards compression therapy      | 24 |
| Dauerbinde® K Long Stretch Bandage             | 22 |
| Durelast® Very Short Stretch Bandage           | 18 |
| Idealbinde™ Short Stretch Bandage              | 17 |
| Komprex® Foam Rubber Pads, Bandage, Sheet      | 14 |
| Komprex® II Channel Foam Sheet                 | 15 |
| Lenkelast® Medium Stretch Bandage              | 19 |
| Mollelast® Conforming Bandage                  | 10 |
| Panelast® Short Stretch Adhesive Bandage       | 29 |
| Porelast® Short Stretch Adhesive Bandage       | 28 |
| Rosidal® CC Cohesive Compression Bandage       | 23 |
| Rosidal® K Short Stretch Bandage               | 16 |
| Rosidal® Lymphset                              | 20 |
| Rosidal® SC Compression Bandage, soft          | 23 |
| Rosidal® soft Foam Padding Bandage             | 13 |
| Rosidal® sys Short Stretch Compression System  | 30 |
| Rosidal® TCS Two Component System              | 23 |
| Silkafix® Adhesive Tape                        | 21 |
| tg® Cotton Stockinette                         | 8  |
| tg® fix Tubular Net Bandage                    | 31 |
| tg® grip Elasticated Tubular Support Bandage   | 26 |
| tg® shape Shaped Tubular Bandage               | 27 |
| tg® soft Tubular Padding Stockinette           | 9  |
| tg® Tubular Stockinette                        | 7  |
| Transelast® Classic Conforming Bandage         | 11 |
| Unna's Boot, Varicex® Zinc Paste Bandage       | 25 |



# Lymphedema Bandaging Materials



## Multi-layer Bandaging Products for Lymphedema from Lohmann & Rauscher

| Products  |  |   |   |   |   |
|---|--|---|---|---|---|
| Skin Protection   | Fingers / Toes   | Paddings  | Focal Compression (optional)  | Bandage   | Tape  |
| tg® Tubular<br>            | Mollelast®<br>          | Cellona® Synthetic<br> | Komprex® Pads<br>          | Rosidal® K Short Stretch<br> | Silkafix® Tape<br> |
| tg® Cotton Stockinette<br> | Transelast® Classic<br> | Rosidal® soft<br>      | Komprex® Rolls & Sheet<br> | Idealbinde - Short<br>       |   |
| tg® soft<br>               |  |   | Komprex® II<br>            | Lenkelast® - Medium<br>      |   |
|   |  |   |   | Durelast® - Extra Short<br>  |   |

### tg® Tubular Bandage, tg® soft or tg® Cotton Stockinette

These Stockinettes are used to protect the skin and absorb excess perspiration and oil making this the first layer applied.

### Mollelast® or Transelast® Classic Conforming Bandage

These conforming bandages (4 or 6 cm widths) are used to wrap individual fingers or toes, which can be folded in half. The wide sizes, 8, 10 or 12 cm widths, can be used to keep foam padding in place before short stretch bandages are applied.

### Cellona® Synthetic Padding or Rosidal® soft Foam Bandage

These padding products are applied just prior to the actual short stretch compression bandages. They cushion the limb and prevent sharp indentations and/or irritations of the skin. Furthermore, they serve to distribute the pressure evenly over the area of an arm or a leg. For maximum cushioning in skin creases or bony prominences (e.g. popliteal fossa, ankle, wrist, etc.) they can be used as a double or triple layer.

### Rosidal® K Short Stretch Compression Bandage

Rosidal K is 100% cotton. Its stretch, which is necessary to achieve required compression, is due solely to its weave. Because of the specialized bandaging technique and pressure required in lymphedema, we offer several widths and lengths.

### Komprex® High Density Foam Rubber

The use of pre-made or custom cut foam pads are recommended for filling anatomical indentations surrounding bony prominences (e.g. ankles, wrist). Pieces of foam rubber are cut to size and placed in the indentation, which results in more even pressure distribution. Foam rubber pads are also used to achieve localized pressure increases in fibrotic areas.

### Komprex® II Foam, wavy

Narrow strips of pre-made foam embedded in skin-friendly, non-woven material are used to create this breathable, non-slip padding. The special wavy surface structure of Komprex® II creates a firm yet gentle micro-massaging effect to help break up fibrotic tissue.

### Silkafix® Adhesive Tape

This tape is coated with an adherent, skin-friendly polycrylate adhesive and can be removed without leaving any residue on the skin or bandages.



# Recommended Materials for Upper Extremity Lymphedema Bandaging

- Low pH lotion (not available through L&R)
- One tubular bandage  
Size: Varies by product selection / patient  
**tg® Tubular Bandage, tg® soft or tg® Cotton Stockinette**
- Two conforming bandages  
Width: 4 or 6 cm (can be folded in half)  
**Mollelast® or Transelast® Classic**
- Two padding bandages  
Width: 10 cm  
**Cellona® Synthetic Padding or Rosidal® soft**
- High density foam rubber  
**Komprex® Foam Rubber**
- Four to five short stretch compression bandages  
Width: one 6 cm, one 8 cm, two or three 10 cm  
**Rosidal® K Short Stretch Bandage**
- One roll of tape  
**Silkafix® Adhesive Tape**



**Note:** Compression therapy is an essential component of lymphedema management. Bandaging materials will vary by patient. This list is only a suggestion to provide basic information and should not replace the expertise of a trained therapist.

# Recommended Materials for Lower Extremity Lymphedema Bandaging

- Low pH lotion (not available through L&R)
- One tubular bandage  
Size varies by product selection  
**tg® Tubular Bandage, tg® soft or tg® Cotton Stockinette**
- One conforming bandage  
Width: 4 or 6 cm (to be folded in half)  
**Mollelast® or Transelast® Classic**
- Two to four padding bandages  
Width: 10 cm and 15 cm  
**Cellona® Synthetic Padding or Rosidal® soft**
- High density foam rubber  
**Komprex® Foam Rubber**
- Nine to twelve short stretch compression bandages  
Width: one 6 cm, one 8 cm, three or four 10 cm and four to six 12 cm  
**Rosidal® K Short Stretch Bandage**
- One roll of tape  
**Silkafix® Adhesive Tape**

If bandaging of the hip is indicated:

- Two to three **Idealbinde™ Short Stretch Bandages**

**Note:** Compression therapy is an essential component of lymphedema management. Bandaging materials will vary by patient. This list is only a suggestion to provide basic information and should not replace the expertise of a trained therapist.



# tg® Tubular Bandage

## Tubular Stockinette

### Field of Application

- First layer next to skin (lymphology and phlebology)
- Enclosing foam padding structures
- Dressing retention
- Over zinc paste bandages
- Used to make homemade chip foam bags

### Features/Benefits

- Close-meshed, seamlessly knitted tubular bandage
- Lightweight and breathable
- Absorbs excess perspiration
- Stretches widthwise approximately four times its original width
- Can be narrowed down again by pulling lengthwise
- Crease-free application; conforms comfortably to body contours
- Not made with natural rubber latex

### Product composition

cotton, viscose



### tg® Tubular Bandage

| Size  | Item No. | Shipping Units |
|---|----------|----------------|
| <b>20 m (22 yd) length, rolled up in a cardboard dispenser</b>            |          |                |
| 1 (1.4 cm) / 0.55" for finger and toe                                     | 24000    | 1/bx, 20 bx/cs |
| 2 (2.3 cm) / 0.91" for fingers, toes, child's hand                        | 24001    | 1/bx, 20 bx/cs |
| 3 (3.0 cm) / 1.20" for several fingers, child's arm                       | 24002    | 1/bx, 10 bx/cs |
| 5 (5.5 cm) / 2.20" for arm, lower leg, child's leg                        | 24003    | 1/bx, 10 bx/cs |
| 6 (6.5 cm) / 2.60" for arm, leg (medium)                                  | 24004    | 1/bx, 10 bx/cs |
| 7 (7.0 cm) / 2.80" for arm (large), leg, child's head                     | 24005    | 1/bx, 18 bx/cs |
| 9 (8.5 cm) / 3.30" for thigh (large), head, armpit                        | 24006    | 1/bx, 12 bx/cs |
| 12 (12.0 cm) / 4.70" for thigh (extra large), head (large), child's trunk | 24007    | 1/bx, 12 bx/cs |
| <b>10 m (11 yd) length, rolled up in a cardboard dispenser</b>            |          |                |
| K1 (16.0 cm) / 6.30" for trunk (small)                                    | 24030    | 1/bx, 10 bx/cs |
| K2 (21.0 cm) / 8.27" for trunk (large)                                    | 24031    | 1/bx, 10 bx/cs |



# tg<sup>®</sup> Cotton Stockinette

## Tubular Stockinette

### Field of Application

- First layer next to skin (lymphology and phlebology)
- Over zinc paste bandages

### Features/Benefits

- 100% cotton, close-weave, seamlessly knitted tubular stockinette
- Very soft, supple, skin-friendly protection
- Excellent for sensitive skin
- Absorbs excess perspiration
- Stretches widthwise approximately twice its original width
- Retains shape well

### Product composition

100% raw cotton, unbleached

- Not made with natural rubber latex



### tg<sup>®</sup> Cotton Stockinette

| Width   | Item No. | Shipping Units    |
|---|----------|-------------------|
| 25 m (27.5 yd) length, individually sealed in a plastic bag |          |                   |
| 6 cm (2.4")   | 14820    | 1/bag, 10 bags/cs |
| 8 cm (3.2")   | 14821    | 1/bag, 10 bags/cs |
| 10 cm (4.0")  | 14822    | 1/bag 5 bags/cs   |
| 12 cm (4.7")  | 14823    | 1/bag, 5 bags/cs  |
| 15 cm (6.0")  | 14824    | 1/bag, 6 bags/cs  |
| 20 cm (7.9")  | 14825    | 1/bag, 3 bags/cs  |
| 25 cm (9.8")  | 14826    | 1/bag, 3 bags/cs  |
| 30 cm (11.8")   | 14827    | 1/bag, 3 bags/cs  |

# tg<sup>®</sup> soft

## Tubular Padding Stockinette

### Field of Application

- First layer next to skin (lymphology and phlebology)
- For light padding and skin protection
- Over zinc paste bandages

### Features/Benefits

- Soft, comfortable terry cloth tubular bandage
- Absorbs high amounts of perspiration and oils
- Excellent for fragile skin
- Provides extra comfort and padding
- Stretches widthwise approximately five times its original width

### Product composition

85% cotton (unbleached), 15% spandex

- Not made with natural rubber latex



### tg<sup>®</sup> soft

Tubular Padding Bandage

|   | Size | Unstretched Measurement       | Stretched Measurement         | Item No. | Shipping Units |
|---|------|-------------------------------|-------------------------------|----------|----------------|
| 10 m (11 yd) length, rolled up in a cardboard dispenser |      |                               |                               |          |                |
| for forearm, child's arm, child's leg                   | S    | 5.5 - 6.5 cm<br>(2.2" - 2.6") | 30 - 32 cm<br>(11.8" - 12.6") | 13200    | 1/bx, 8 bx/cs  |
| for lower leg, arm                                      | M    | 8 - 10 cm<br>(3.2" - 4.0")    | 41 - 48 cm<br>(16.1" - 19.0") | 13201    | 1/bx, 6 bx/cs  |
| for leg   | L    | 13 - 15 cm<br>(5.1" - 5.9")   | 56 - 66 cm<br>(22.1" - 26.0") | 13202    | 1/bx, 4 bx/cs  |

# Mollelast®

## Conforming Bandage

### Field of Application

- Compression bandaging for lymphedema
- Fingers, Toes and Genital lymphedema
- Fixing foam padding
- Dressing retention

### Features/Benefits

- Comfortable, open-weave gauze bandage
- Ideal for bandaging fingers and toes
- Conforms well to body contours
- Air permeable
- Wider sizes are excellent to secure foam padding from moving under bandages
- Extensibility is approximately 100%

### Product composition

56% viscose, 44% polyamide

- Not made with natural rubber latex



### Mollelast®

#### Conforming Bandage

| Width                         | Item No. | Shipping Units        |
|-------------------------------|----------|-----------------------|
| 4 m (4.4 yd) stretched length |          |                       |
| 4 cm (1.6")                   | 19410    | 20 rolls/bx, 26 bx/cs |
| 6 cm (2.4")                   | 19411    | 20 rolls/bx, 18 bx/cs |
| 8 cm (3.2")                   | 19412    | 20 rolls/bx, 25 bx/cs |
| 10 cm (4.0")                  | 19413    | 20 rolls/bx, 18 bx/cs |
| 12 cm (4.7")                  | 19414    | 20 rolls/bx, 16 bx/cs |

### Mollelast®





# Transelast® Classic

## Conforming Bandage

### Field of Application

- Compression bandaging for lymphedema
- Fingers, Toes and Genital lymphedema
- Dressing retention

### Features/Benefits

- Soft, comfortable, close-weave gauze bandage
- Ideal for bandaging fingers and toes
- Conforms well to body contours
- Air permeable
- Extensibility is approximately 130%

### Product composition

60% polyamide, 40% viscose

- Not made with natural rubber latex



### Transelast® Classic

#### Conforming Bandage

| Width                                | Item No. | Shipping Units        |
|--------------------------------------|----------|-----------------------|
| <b>4 m (4.4 yd) stretched length</b> |          |                       |
| 4 cm (1.6")                          | 90400    | 20 rolls/bx, 33 bx/cs |
| 6 cm (2.4")                          | 90401    | 20 rolls/bx, 21bx/cs  |

### Transelast® Classic



# Cellona®

## Synthetic Padding

### Field of Application

- Padding under compression bandaging (lymphology and phlebology)

### Features/Benefits

- Cushions the limb
- Prevents sharp indentations and/or irritations of the skin
- Ensures even distribution of pressure from compression bandages
- Material's self-adhesion facilitates application and helps prevent slippage
- Soft and elastic; provides comfort and excellent conformability around body contours
- Helps form a smooth surface
- Air permeable
- Retains fullness with long-term wear

### Box packaging:

Each of the four sizes of padding are packaged in the same size box.

- Same size boxes are easy to stack and store
- Box labeling easily identifies product size
- Professional looking, self contained box instead of loose rolls

### Product composition

100% polyester fibers

- Not made with natural rubber latex



### Cellona®

#### Synthetic Padding

| Size  | Item No. | Shipping Units       |
|---|----------|----------------------|
| individually wrapped in paper sleeve inside box |          |                      |
| 4 cm x 3 m (1.6" x 3.3 yd)                      | 10682    | 9 rolls/bx, 10 bx/cs |
| 6 cm x 3 m (1.6" x 3.3 yd)                      | 10683    | 6 rolls/bx, 10 bx/cs |
| 10 cm x 3 m (4" x 3.3 yd)                       | 10684    | 4 rolls/bx, 10 bx/cs |
| 15 cm x 3 m (6" x 3.3 yd)                       | 10685    | 3 rolls/bx, 10 bx/cs |

### Cellona®

#### Synthetic Padding

| Size   | Item No. | Shipping Units |
|--|----------|----------------|
| individually wrapped in paper sleeve inside bulk box |          |                |
| 4 cm x 3 m (1.6" x 3.3 yd)                           | 10690    | 108 rolls/bx   |
| 6 cm x 3 m (1.6" x 3.3 yd)                           | 10692    | 72 rolls/bx    |
| 10 cm x 3 m (4" x 3.3 yd)                            | 10694    | 48 rolls/bx    |
| 15 cm x 3 m (6" x 3.3 yd)                            | 10695    | 36 rolls/bx    |

# Rosidal® soft

## Foam Padding Bandage

### Field of Application

- Padding under compression bandaging (lymphology and phlebology)

### Features/Benefits

- Reusable; washable up to 50 times
- Ensures even distribution of pressure from compression bandages
- Interlocking open pores make compression bandages slip resistant
- Durable and tear resistant
- Provides some compression and shaping of the limb when pulled slightly
- Highly breathable and permeable to moisture
- Washing temperature up to 95°C/203°F (temperatures from 104°-140°F/ 40°-60°C are recommended)
- Maintains elasticity and cushioning effect even after washing

### Product composition

100% polyurethane foam

- Product and clips are not made with natural rubber latex



### Rosidal® soft

Foam Padding Bandage

| Size  | Item No. | Shipping Units |
|---|----------|----------------|
| individually packed in cardboard box              |          |                |
| 10 cm x 0.3 cm x 2.5 m<br>(4" x 0.12" x 2.7 yd)   | 23110    | 1/bx, 30 bx/cs |
| 10 cm x 0.4 cm x 2.5 m<br>(4" x 0.16" x 2.7 yd)   | 23111    | 1/bx, 24 bx/cs |
| 12 cm x 0.4 cm x 2.5 m<br>(4.7" x 0.16" x 2.7 yd) | 23112    | 1/bx, 16 bx/cs |
| 15 cm x 0.4 cm x 2.5 m<br>(6" x 0.16" x 2.7 yd)   | 23113    | 1/bx, 14 bx/cs |



# Kompres<sup>®</sup>

## Foam Rubber Pads, Bandage, Sheet

### Field of Application

- Compression in treatment of edema (lymphology and phlebology)
- Increase local pressure (in combination with other compression products)
- Soften fibrosis
- Fill anatomical depressions

### Features/Benefits

- Reusable, closed-cell, dense foam
- Easily cuts to required shape and size
- Available in ready-to-use shapes, rolls and sheet
- Permeable to air and water vapor

### Product composition

100% synthetic latex

- Not made with natural rubber latex

**Note:** If Kompres<sup>®</sup> Foam is cut to a specific shape, the edges should be chamfered to prevent pressure spots.



### Kompres<sup>®</sup> Foam Rubber Pads

individually packed in a cardboard box

| Size            | Thickness    | Size                      | Item No. | Shipping Units |
|-----------------|--------------|---------------------------|----------|----------------|
| 00 oval         | 10 mm (0.4") | 9 x 6.5 cm (3.5" x 2.6")  | 22300    | 1/bx, 60 bx/cs |
| 0 kidney-shaped | 10 mm (0.4") | 9 x 5.0 cm (3.5" x 2.0")  | 22301    | 1/bx, 75 bx/cs |
| 1 kidney-shaped | 12 mm (0.5") | 12 x 6.5 cm (4.8" x 2.6") | 22302    | 1/bx, 30 bx/cs |

### Kompres<sup>®</sup> Foam Rubber Bandage

rolled, individually packed in a cardboard box

| Thickness    | Size                       | Item No. | Shipping Units |
|--------------|----------------------------|----------|----------------|
| 5 mm (0.2")  | 8 cm x 2 m (3.2" x 2.2 yd) | 22313    | 1/bx, 30 bx/cs |
| 10 mm (0.4") | 8 cm x 2 m (3.2" x 2.2 yd) | 22314    | 1/bx, 15 bx/cs |

### Kompres<sup>®</sup> Foam Rubber Sheet

rolled, in plastic bag

| Thickness    | Size                         | Item No. | Shipping Units    |
|--------------|------------------------------|----------|-------------------|
| 10 mm (0.4") | 100 cm x 50 cm (39.4" x 20") | 22320    | 1/bag, 10 bags/cs |

# Komprex® II

## Channel Foam Sheet

### Field of Application

- Compression in treatment of edema (lymphology and phlebology)
- Increase local pressure (in combination with other compression products)
- Soften fibrosis

### Features/Benefits

- Special surface creates a firm, yet gentle micro-massaging effect
- Narrow strips of open cell foam material are surrounded by a skin-friendly, non-woven, breathable, non-slip padding
- Provides channeling effect
  - Loosens hardened and fibrotic tissue
  - Better air flow for greater patient comfort and compliance
- Washable (up to 30°C/86°F) and reusable.

### Product composition

Narrow strips of polyurethane foam rubber embedded in a polyester, non-woven fabric

- Not made with natural rubber latex

**Note:** To avoid obstructions of lymphatic flow, the foam line must always be applied in a direction lengthwise to the extremity with the soft, comfortable side of the material facing the body.



### Komprex® II

#### Foam Sheet

| Size   | Item No. | Shipping Units   |
|--|----------|------------------|
| individually wrapped in paper and a resealable plastic bag |          |                  |
| 65 cm x 65 cm (26" x 26")                                  | 16985    | 1/bag, 6 bags/cs |

# Rosidal® K

## Short Stretch Bandage

### Field of Application

- Strong compression bandaging (lymphology and phlebology)
- Support and stress relief in acute care and sports medicine

### Features/Benefits

- 100% cotton
- Reusable; washable up to 50 times
- Provides high working pressure and low resting pressure well tolerated even at night
- Well tolerated, even at night
- Bandage can be worn at rest and during exercise
- Distinct textile structure helps prevent bandage layers from slipping
- Very skin-friendly
- Air permeable
- Simplify application of thigh/large leg by using the double (10 m) length bandage
- Extensibility is approximately 90%
- Washing temperature up to 95°C/203°F

### Product composition

100% cotton

- Product and clips are not made with natural rubber latex

### Rosidal® K



### Rosidal® K

Short Stretch Bandage

| Width   | Item No. | Shipping Units  |
|---|----------|-----------------|
| 5 m (5.5 yd) length, individually packed in a cardboard box |          |                 |
| 4 cm (1.6")   | 90684    | 1/bx, 20 bx/cs  |
|   | 22199    | 1/bx, 125 bx/cs |
| 6 cm (2.4")   | 90685    | 1/bx, 20 bx/cs  |
|   | 22200    | 1/bx, 108 bx/cs |
| 8 cm (3.2")   | 90686    | 1/bx, 20 bx/cs  |
|   | 22201    | 1/bx, 108 bx/cs |
| 10 cm (4.0")  | 90687    | 1/bx, 20 bx/cs  |
|   | 22202    | 1/bx, 90 bx/cs  |
| 12 cm (4.7")  | 90688    | 1/bx, 20 bx/cs  |
|   | 22203    | 1/bx, 54 bx/cs  |
| 10 m (11 yd) length, individually packed in a cardboard box |          |                 |
| 10 cm (4.0")  | 90689    | 1/bx, 20 bx/cs  |
|   | 22206    | 1/bx, 60 bx/cs  |
| 12 cm (4.7")  | 90690    | 1/bx, 20 bx/cs  |
|   | 22205    | 1/bx, 40 bx/cs  |



# Idealbinde™

## Short Stretch Bandage

### Field of Application

- Strong compression bandaging (lymphology and phlebology)
- Edema in the chest and abdomen
- Bandage stabilization on the thigh and abdomen
- Fixing foam padding applications
- Support and stress relief in cases of distortion or contusion

### Features/Benefits

- 100% cotton, wide-width bandage
- Firm grip fabric structure ensures bandage layers do not slip
- Skin-friendly
- Air permeable
- Extensibility is approximately 90%
- Washable up to 95°C/203°F

### Product composition

100% cotton

- Product and clips are not made with natural rubber latex



### Idealbinde™

#### Short Stretch Bandage

| Width   | Item No. | Shipping Units       |
|---|----------|----------------------|
| 5 m (5.5 yd) stretched length, individually wrapped |          |                      |
| 15 cm (6.0")  | 14105    | 10 rolls/bx, 6 bx/cs |
| 20 cm (7.9")  | 14106    | 35 rolls/cs          |

# Durelast®

## Very Short Stretch Bandage

### Field of Application

- Very strong compression bandaging (lymphology and phlebology)
- Support and stress relief in acute care and sports medicine

### Features/Benefits

- Firm grip fabric texture ensures bandage layers hold securely
- Very skin-friendly
- Air permeable
- Extensibility is approximately 45%
- Washable up to 95°C/203°F

### Product composition

66% cotton, 34% polyamide

- Product and clips are not made with natural rubber latex



### Durelast®

#### Very Short Stretch Bandage

| Width   | Item No. | Shipping Units  |
|---|----------|-----------------|
| 5 m (5.5 yd) stretched length, individually packed in a cardboard box |          |                 |
| 6 cm (2.4")   | 22100    | 1/bx, 108 bx/cs |
| 8 cm (3.2")   | 22101    | 1/bx, 108 bx/cs |
| 10 cm (4.0")  | 22102    | 1/bx, 90 bx/cs  |
| 12 cm (4.7")  | 22103    | 1/bx, 54 bx/cs  |

Durelast®



# Lenkelast®

## Medium Stretch Bandage

### Field of Application

- Compression bandaging (lymphology and phlebology)
- Genital lymphedema
- Pediatric lymphedema
- Breast bandaging
- Fixing layer over underlying foam padding
- Dressing retention
- Immobilization of limbs

### Features/Benefits

- Soft, comfortable, easy to use all purpose bandage
- Conforms well to all body contours
- Very skin-friendly and breathable
- Durable elasticity for a wide variety of applications
- Extensibility is approximately 130%

### Product composition

64% cotton, 35% polyamide,  
1% spandex

- Not made with natural rubber latex



### Lenkelast®

#### Medium Stretch Bandage

| Width   | Item No. | Shipping Units        |
|---|----------|-----------------------|
| 5 m (5.5 yd) stretched length, individually wrapped |          |                       |
| 6 cm (2.4")   | 19460    | 10 rolls/bx, 25 bx/cs |
| 8 cm (3.2")   | 19461    | 10 rolls/bx, 20 bx/cs |
| 10 cm (4.0")  | 19462    | 10 rolls/bx, 24 bx/cs |
| 12 cm (4.7")  | 19463    | 10 rolls/bx, 20 bx/cs |
| 15 cm (6.0")  | 19464    | 10 rolls/bx, 16 bx/cs |
| 20 cm (7.9")  | 19465    | 10 rolls/bx, 10 bx/cs |

# Rosidal® Lymphset

## Complete kit for Lymphedema treatment

### Field of Application

- Professional lymphological compression bandaging of legs and arms

### Features/Benefits

- Convenient, prepackaged kit saves time and simplifies the order process and inventory management
- Contains all the necessary component required for a professional multi-layer compression bandage
- Kits contain two complete sets of bandaging material - one to wear and one to wash for uninterrupted usage
- Requires no preparation time for specific patient care. It's ready for use at any time.



Item No. 49103 pictured

| Components  | Deluxe Double Arm<br>Item No. 49101   | Deluxe Double Leg<br>Item No. 49103  |
|---|---|--|
| <b>tg® Tubular Bandage</b><br>For skin protection under decongestive compression bandage  | Size 6<br>(1x)  | Size 9<br>(1x)   |
| <b>Mollelast® Conforming Bandage</b><br>Comfortably soft, highly elastic conforming bandage for the compression of fingers and toes in lymphedema bandaging                       | 6 cm x 4 m (2.4" x 4.4 yd)<br>(20x)   | 4 cm x 4 m (1.6" x 4.4 yd)<br>(20x)  |
| <b>Cellona® Synthetic Padding</b><br>Padding under decongestive compression bandages for even distribution of pressure and for the prevention of constrictions and pressure spots | 10 cm x 3 m<br>(4" x 3.3 yd)<br>(4x)  | 10 cm x 3 m (4" x 3.3 yd)<br>(4x)<br>15 cm x 3 m (6" x 3.3 yd)<br>(6x)   |
| <b>Rosidal® K Short Stretch Bandage</b><br>100% cotton, for strong compression of extremities in lymphedema treatment   | 6 cm x 5 m (2.4" x 5.5 yd)<br>(2x)<br>8 cm x 5 m (3.2" x 5.5 yd)<br>(2x)<br>10 cm x 5 m (4" x 5.5 yd)<br>(4x) | 6 cm x 5 m (2.4" x 5.5 yd)<br>(2x)<br>8 cm x 5 m (3.2" x 5.5 yd)<br>(2x)<br>10 cm x 5 m (4" x 5.5 yd)<br>(4x)<br>12 cm x 5 m (4.7" x 5.5 yd)<br>(4x) |
| <b>Idealbinde Short Stretch Bandage</b><br>For strong compression and bandage stabilization of the thigh and abdomen  |   | 15 cm x 5 m (6" x 5.5 yd)<br>(2x)  |
| <b>Komprex® Foam Rubber Pad</b><br>For padding and more even pressure distribution in the ankle joint   |   | Kidney shaped<br>9 cm x 5 cm (3.5" x 2")<br>(2x)   |
| <b>Silkafix® Adhesive Tape</b><br>For securing lymphedema bandages  | 1.25 cm x 5 m (½" x 5.5 yd)<br>(1x)   | 1.25 cm x 5 m (½" x 5.5 yd)<br>(1x)  |



# Silkafix®

## Adhesive Tape

### Field of Application

- Fastens bandages, dressings, catheters, cannulas, drains or probes

### Features/Benefits

- Strong, durable, satiny tape adheres well, yet removes easily
- High performance adhesive leaves no residue on skin or bandages
- Easily tears lengthwise and width-wise for customized application; no scissors needed
- Particularly suitable for sensitive skin
- Skin-friendly, hypoallergenic
- White artificial silk, non-elastic tape
- Water and dirt repellent
- Quick and easy application; dispenses from a quality, protected, plastic spool

### Product composition

rayon taffeta, white, 100% acetate; skin friendly polyacrylate adhesive

- Not made with natural rubber latex



**Silkafix®**  
Adhesive Tape

| Width  | Item No. | Shipping Units |
|--|----------|----------------|
| 5 m (5.5 yd) roll on a plastic dispensing reel |          |                |
| 1.25 cm (0.5")                                 | 34324    | 18 rolls/bx    |
| 2.5 cm (1.0")                                  | 34325    | 12 rolls/bx    |
| 5 cm (2.0")                                    | 34326    | 6 rolls/bx     |



# Bandage Roller

- Quick bandage re-rolling with manually operated device
- Bandage is placed into the slot and draped over the tension bar
- Easy to use; free hand guides bandage while controlling tension
- Roller can be mounted to any flat surface with enclosed bracket
- After usage, bandage roller can be taken out of the bracket holder and stored

## Bandage Roller

| Item No.                               | Shipping Units |
|--|----------------|
| individually packaged in a plastic bag |                |
| CR3000                                 | 1 piece        |

# Dauerbinde® K

## Long Stretch Bandage

### Field of Application

- Compression bandaging (phlebology)
- Support and stress relief to injured tendon, ligament or muscle
- Dressing retention
- Limb immobilization

### Features/Benefits

- Easy to apply and conforms well to body contours without creases or folds
- Compressive power can be easily adjusted due to long stretch properties
- Durable textured fabric structure
- Helps stabilize the compression bandage and keeps it from slipping when applied as a final layer
- Extensibility is approximately 180%
- Washable up to 30°C/86°F

### Product composition

86% cotton, 7% polyamide, 7% spandex

- Product and clips are not made with natural rubber latex

Note: Due to the high resting pressure level, the bandage should not be worn overnight or during extended periods of rest.



### Dauerbinde® K



### Dauerbinde® K

#### Long Stretch Bandage

| Width   | Item No. | Shipping Units  |
|---|----------|-----------------|
| 7 m (7.7 yd) stretched length, individually packed in a cardboard box |          |                 |
| 6 cm (2.4")   | 22000    | 1/bx, 100 bx/cs |
| 8 cm (3.2")   | 22001    | 1/bx, 100 bx/cs |
| 10 cm (4.0")  | 22002    | 1/bx, 100 bx/cs |
| 12 cm (4.7")  | 22003    | 1/bx, 50 bx/cs  |
| 20 cm (7.9")  | 22004    | 1/bx, 20 bx/cs  |

# Rosidal® TCS

## Two Component System

### Field of Application

- Venous leg ulcer
- Edema
- Concomitant conditions

ABPI above 0.8 or with an ABPI of 0.5-0.8 only under medical control

### Features/Benefits

- The two components of Rosidal TCS work together to deliver effective levels of working compression and low levels of resting compression
- With the Safe-Loc system, both components are applied at full stretch: pull to 100% stretch
- Comfort and safety built in with the Safe-Loc System
- First soft compression bandage layer combines padding and skin protection in a single bandage
  - Acts as a cushion to evenly distribute pressure
  - Integrated fabric lining protects sensitive skin from irritation
- Durable short stretch cohesive bandage layer
  - Adhesive helps secure the compression bandage
- System is thin enough to allow the use of normal footwear to help encourage compliance



### Rosidal® TCS

#### Two Component System

| System                                     | Size                         | Item No. | Shipping Units |
|--|------------------------------|----------|----------------|
| packed together in a cardboard box         |                              |          |                |
| Rosidal SC<br>compression bandage, soft    | 4" x 3.8 yd. (10 cm x 3.5 m) | 26484    | 1/bx, 14 bx/cs |
| Rosidal CC<br>cohesive compression bandage | 4" x 6.6 yd. (10 cm x 6.0 m) |          |                |

### Product composition

#### Rosidal SC

Polyurethane foam, with woven textile material made from 60% polyamide and 40% cotton

- Not made with natural rubber latex

#### Rosidal CC

51% cotton, 39% polyamide, 10% polyurethane, with a hypo-allergenic synthetic microcoating

- Not made with natural rubber latex

### Rosidal® SC and Rosidal CC

| Individual Components                      | Size                         | Item No. | Shipping Units |
|--|------------------------------|----------|----------------|
| individually packed in a cardboard box     |                              |          |                |
| Rosidal SC<br>compression bandage, soft    | 4" x 3.8 yd. (10 cm x 3.5 m) | 33503    | 1/bx, 24 bx/cs |
| Rosidal CC<br>cohesive compression bandage | 4" x 6.6 yd. (10 cm x 6.0 m) | 33499    | 1/bx, 24 bx/cs |



# Compression Bandaging for Lymphedema Treatment


by Guenter Klose

- Provides a complete overview of lymphedema, the role and function of bandaging, and specific techniques to benefit the lymphedema patient

Revision date 2014

## Compression Bandaging for Lymphedema Treatment

| Item No.      | Shipping Units |
|---------------|----------------|
| 7894603/R0714 | 1 piece        |



# Contributions towards compression therapy

by Prof. Dr. H. Partsch

- Provides an outstanding overview of venous disorders and the use of compression bandaging in their treatment

Revision date 2014

## Contributions towards compression therapy

| Item No. | Shipping Units |
|----------|----------------|
| SL1035   | 1 piece        |

# Unna's Boot, Varicex® Zinc Paste Bandage

## Field of Application

- Semi-rigid treatment of venous insufficiency, leg ulcers, edemas, thrombophlebitis
- Secondary treatment of fractures and dislocation

## Features/Benefits

- Ready-to-use zinc oxide impregnated bandages
- Easy application with easy unrolling; firm center core prevents roll from collapsing thus minimizing constrictions
- Rolls are individually wrapped in a foil pouch to keep bandage fresh
- Two varieties available
  - **Varicex® S** - easily conforms to body contours; knitted bandage with widthwise and lengthwise stretch
  - **Varicex® F** - gauze bandage with low moisture content for rapid drying

## Product composition

**Varicex® S** - viscose/polyamide bandage impregnated with zinc paste

- Not made with natural rubber latex

**Varicex® F** - cotton and viscose ribbon gauze impregnated with zinc paste

- Not made with natural rubber latex



## Varicex® S Unna's Boot

| Width   | Item No. | Shipping Units         |
|---|----------|------------------------|
| 10 m (11 yd.) stretched length, individually sealed |          |                        |
| 10 cm (4")  | 20632    | 10 pouches/bx, 3 bx/cs |

## Varicex® F Unna's Boot

| Width   | Item No. | Shipping Units         |
|---|----------|------------------------|
| 10 m (11 yd.) unstretched length, individually sealed |          |                        |
| 10 cm (4")  | 20602    | 10 pouches/bx, 5 bx/cs |

# tg® grip

## Elasticated Tubular Support Bandage

### Field of Application

- Light compression for effective management of strains, sprains, and soft tissue injuries
- Support and relief of ligaments, tendons and joints
- Dressing retention

### Features/Benefits

- Seamlessly knitted tube adjusts to body contours and evenly distributes pressure
- Comfortable, skin-friendly, breathable elasticated cotton construction
- Quick and easy application
- Remove and reapply without clips, pins or tapes
- Widest range of sizes available for most applications
- Washable (60°C/140°F) and reusable

### Product composition

cotton, elastodiene (natural rubber latex) and polyamide

**Tip:** Compression can be enhanced with a second layer of tg® grip or by using a smaller size of tg® grip.

tg® grip



tg® grip  
Elasticated Tubular Support Bandage

|   | Size  | Width         | Item No. | Shipping Units |
|---|---|---------------|----------|----------------|
| 10 m (11 yd.) roll in a cardboard dispenser |   |               |          |                |
| A   | Children's hands, children's forearms                               | 1¾" (4.5 cm)  | 24120    | 1/bx, 24 bx/cs |
| B   | Small hand/wrist  | 2⅓" (6.0 cm)  | 24121    | 1/bx, 16 bx/cs |
| C   | Medium hand/wrist<br>small elbow/ankle                              | 2¾" (6.5 cm)  | 24122    | 1/bx, 16 bx/cs |
| D   | Large hand/wrist, medium elbow/ankle, small leg/knee                | 3" (7.5 cm)   | 24123    | 1/bx, 16 bx/cs |
| E   | X-Large hand/wrist, large elbow/ankle, medium leg/knee, small thigh | 3½" (8.75 cm) | 24124    | 1/bx, 12 bx/cs |
| F   | X-Large elbow/ankle, large leg/knee, medium thigh                   | 4" (10 cm)    | 24125    | 1/bx, 12 bx/cs |
| G   | X-Large leg/knee, large thigh                                       | 4¾" (12 cm)   | 24126    | 1/bx, 8 bx/cs  |
| J   | X-Large thigh, small trunk  | 7" (17.5 cm)  | 24127    | 1/bx, 6 bx/cs  |
| K   | Medium trunk  | 8½" (21.5 cm) | 24128    | 1/bx, 4 bx/cs  |
| L   | Large trunk   | 13" (32.5 cm) | 24129    | 1/bx, 4 bx/cs  |
| M   | X-Large trunk   | 15" (37.5 cm) | 24130    | 1/bx, 2 bx/cs  |



# tg<sup>®</sup> shape

## Shaped Tubular Bandage

### Field of Application

- Light compression and support

### Features/Benefits

- Anatomically shaped tubular bandage to provide compression and fixation along the leg
- Easy to apply and reapply without the use of clips, pins or tape
- Skin friendly, breathable elastic cotton construction
- Washable (30°C/86°F) and reusable
- Allows freedom of movement
- Quick and easy size selection
- Two convenient styles available:
  - Below knee
  - Full Leg
- Seamless Knit Design

### Product composition

cotton, elastodiene (natural rubber latex) and polyamide



### tg<sup>®</sup> shape

#### Shaped Tubular Bandage

| Size                 | Calf Circumference      | Item No. | Shipping Units |
|----------------------|-------------------------|----------|----------------|
| individually wrapped |                         |          |                |
| Full leg S           | 12½ - 14¼" (32 - 36 cm) | 88900    | 10/bx          |
| Full leg M           | 13¾ - 15¼" (35 - 39 cm) | 88901    | 10/bx          |
| Full leg L           | 15 - 16½" (38 - 42 cm)  | 88902    | 10/bx          |
| Full leg XL          | 16¼ - 17¾" (41 - 45 cm) | 88903    | 10/bx          |
| Below knee S         | 12½ - 14¼" (32 - 36 cm) | 88904    | 10/bx          |
| Below knee M         | 13¾ - 15¼" (35 - 39 cm) | 88905    | 10/bx          |
| Below knee L         | 15 - 16½" (38 - 42 cm)  | 88906    | 10/bx          |

### tg<sup>®</sup> shape



# Porelast®

## Short Stretch Adhesive Bandage

### Field of Application

- Compression bandaging in phlebology
- Functional bandaging in sports medicine
- Support and relief with distortions, contusions and dislocations

### Features/Benefits

- Lengthwise stretchable woven, 100% cotton bandage
- Fabric with porously applied adhesive is permeable to air and water vapor
- X-ray transparent
- Extensibility approximately 80%

### Product composition

100% cotton with synthetic rubber adhesive (contains zinc oxide)

- Not made with natural rubber latex



### Porelast®

#### Short Stretch Adhesive Bandage

| Width   | Item No. | Shipping Units |
|---|----------|----------------|
| 2.5 m (2.7 yd) unstretched length, individually packed in a cardboard box |          |                |
| 6 cm (2.4")   | 30210    | 1/bx, 50 bx/cs |
| 8 cm (3.2")   | 30211    | 1/bx, 50 bx/cs |
| 10 cm (4.0")  | 30212    | 1/bx, 50 bx/cs |

### Porelast®



# Panelast®

## Short Stretch Adhesive Bandage

### Field of Application

- Compression bandaging in phlebology
- Base layer in tape applications in sports medicine
- Support and relief with distortions, contusions and dislocations
- Follow-up treatment of fractures

### Features/Benefits

- Lengthwise and widthwise stretchable woven, 100% cotton bandage
- Excellent conformability around body contours
- Fabric with porously applied adhesive is permeable to air and water vapor
- X-ray transparent
- Extensibility widthwise approximately 40%
- Extensibility lengthwise approximately 80%

### Product composition

100% cotton with synthetic rubber adhesive (contains zinc oxide)

- Not made with natural rubber latex



### Panelast®

#### Short Stretch Adhesive Bandage

| Width  | Item No. | Shipping Units |
|--|----------|----------------|
| 2.5 m (2.7 yd.) unstretched length, individually packed in a cardboard box |          |                |
| 6 cm (2.4")  | 30270    | 1/bx, 50 bx/cs |
| 8 cm (3.2")  | 30271    | 1/bx, 50 bx/cs |
| 10 cm (4.0")   | 30272    | 1/bx, 50 bx/cs |

### Panelast®





# Rosidal® Sys

## Short Stretch Compression System

**Field of Application**

- Compression therapy of venous leg ulcers

**Features/Benefits**

- Ready-to-use products combined in one handy set
- Cost effective with reusable padding (Rosidal® soft) and compression bandages (Rosidal® K)
- Rosidal® K short stretch bandage provides high working pressure and low resting pressure making the system tolerable at night
- System can be worn continuously for up to 7 days
- Set saves handling time and eliminates administration of individual components
- Rosidal® soft and Rosidal® K are washable up to 50 times at up to 95°C/203°F



Mollelast® haft contains natural rubber latex

Note: For wound treatment, a suitable dressing can be easily added so the wound is optimally cared for.

### Rosidal® Sys

#### Short Stretch Compression System

| Set         | Item No. | Shipping Units   |
|-------------|----------|------------------|
| Rosidal Sys | 16984    | 1 set, 4 sets/cs |

| Set Components  | Quantity | Size  |
|---|----------|---|
| <b>tg® Tubular Bandage</b><br>For skin protection under the decongestive compression bandage  | 1x       | Size 7  |
| <b>Rosidal® soft Foam Padding Bandage</b><br>Reusable, skin-friendly, latex free foam padding for the even distribution of pressure underneath compression bandage; washable up to 95°C/203°F | 4x       | 10 cm x 2 m x 0.2 cm<br>(4" x 2.2 yd x 0.08")                 |
| <b>Rosidal® K Short Stretch Bandage</b><br>Reusable, 100% cotton short stretch bandage for strong compression of the extremity; washable up to 95°C/203°F                                     | 2x<br>2x | 8 cm x 5 m<br>(3.2" x 5.5 yd)<br>10 cm x 5 m<br>(4" x 5.5 yd) |
| <b>Porofix® Adhesive Tape</b><br>Skin-colored adhesive tape on a roll; coated with synthetic rubber adhesive  | 1x       | 2.5 cm x 5 m<br>(1" x 5.5 yd)                                 |
| <b>Mollelast® haft Cohesive Bandage</b><br>Soft cohesive conforming bandage, lengthwise and widthwise elastic   | 1x       | 8 cm x 20 m<br>(3.2" x 22 yd)                                 |



# tg<sup>®</sup> fix

## Tubular Net Bandage

### Field of Application

- Retention of all types of wound dressings

### Features/Benefits

- Wide-mesh, highly elastic in both directions tubular net
- Fast and easy dressing inspection
- Material is ladderproof; it can be cut open at any point without the cut tearing
- Conformability can be further increased by placing incisions in the right places without incurring the risk of fraying
- Excellent conformability around body contours
- Soft, comfortable and does not restrict movement
- Washable (up to 60°C/140°F) and reusable



### Product composition

71% polyamide,  
29% elastodiene (natural  
rubber latex)

### tg<sup>®</sup> fix Tubular Net Bandage

| Description  | Size | Width<br>Unstretched | Width<br>Stretched | Item No. | Shipping<br>Units |
|--|------|----------------------|--------------------|----------|-------------------|
| 2.5 m (2.7 yd.) stretched length, individually packed in a cardboard box |      |                      |                    |          |                   |
| for finger   | A    | 1.4 cm (0.5")        | 17.5 cm (6.9")     | 24250    | 1/bx, 66 bx/cs    |
| for small extremities<br>(several fingers together/hand/foot)            | B    | 2.7 cm (1.0")        | 40 cm (15.7")      | 24251    | 1/bx, 20 bx/cs    |
| for large extremities  | C    | 4.2 cm (1.7")        | 65 cm (25.6")      | 24252    | 1/bx, 16 bx/cs    |
| for large head, small trunk  | D    | 7.0 cm (3.0")        | 105 cm (41.3")     | 24253    | 1/bx, 16 bx/cs    |
| for large trunk, hip, armpit   | E    | 9.1 cm (3.6")        | 140 cm (55.1")     | 24254    | 1/bx, 12 bx/cs    |

### tg<sup>®</sup> fix





**For more information or to find a dealer near you, contact:**

**L&R USA, Inc.**

3880 W Wheelhouse Rd.

Milwaukee, WI 53208

Phone: 800-279-3863

Fax: 800-289-1903

E-mail: [inquiries@us.LRmed.com](mailto:inquiries@us.LRmed.com)

Website: [www.lohmann-rauscher.us](http://www.lohmann-rauscher.us)

Registered trademarks of Lohmann & Rauscher: Rosidal®, tg®, Komplex®,  
Mollelast®, Transelast®, Cellona®, Durelast®, Lenkelast®, Dauerbinde®,  
Silkafix®, Varicex®, Porelast®, Panelast®