

sigvaris



Compreflex®
INELASTIC COMPRESSION WRAPS

SIGVARIS
GROUP

Simple Solutions for
Venous Disease, Lymphedema & Wound Care

Compreflex Standard

CALF & FOOT

20-50mmHg



When to Use

- Patient requires a calf and foot garment
- Varying ranges of adjustable resting compression is desired (20-50mmHg)
- Non-attached foot piece is desired
- No compression in the liner is desired

Indications

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lipedema



Compreboot® Standard

Variable compression
(see page 7 for details)



Pair of Basic Liners

Non-Compressive
(see page 6 for details)



Black



Beige



Accutab® Devices

See page 4 for details

Compreflex Standard

CALF

20-50mmHg



When to Use

- Patient requires compression in the calf only or a different foot option
- Varying ranges of adjustable resting compression is desired (20-50mmHg)
- No compression in the liner is desired
- No Compreboot needed
- Ability to wear standard footwear over liner is desired

Indications

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lipedema



Pair of Basic Liners

Non-Compressive
(see page 6 for details)



Black



Beige



Accutab Devices

See page 4 for details



Compreflex Transition

CALF

20–50mmHg



When to Use

- Patient requires a compressive liner; the Transition Liner provides 10-15mmHg in the foot and ankle
- Varying ranges of adjustable resting compression is desired (20–50 mmHg)
- Ability to wear standard footwear over liner is desired

Indications

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lipedema



Pair of Transition Liners

10–15 mmHg of mild compression in the foot and ankle.
(see page 6 for details)



Billable under Medicare code A6545-AW



Black



Beige



Accutab Devices

See page 4 for details

Compreflex Complete

CALF

20–50mmHg



When to Use

- Patient requires a compressive liner; the Complete Liner provides 10-15mmHg in the foot, ankle, and calf
- Varying ranges of adjustable resting compression is desired (20–50 mmHg)
- Ability to wear standard footwear over liner is desired

Indications

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lipedema



Pair of Complete Liners

10–15 mmHg of mild compression in the foot, ankle, and calf.
(see page 6 for details)



Billable under Medicare code A6545-AW



Black



Beige



Accutab Devices




See page 4 for details

Easily And Accurately Set Compression With Accutab




Easy To Use

1. Gently lay the Compriflex straps in place without stretching the fabric
2. Attach an Accutab at the end of each strap
3. Pull each strap to the prescribed dosage. Compression ranges are indicated by the lines on the Accutab; 20-30, 30-40, or 40-50 mmHg.



Proven Accuracy For Resting Pressure

Accutab Results: 20-30mmHg



Individual Tests


* Accutab has been rigorously tested according to established guidelines in the International Compression Club Consensus documents^{1,2}. Testing was conducted using PicoPress, which is recognized for its high accuracy and reproducibility³.

Effective Working Pressure

- Inelastic compression reduces leg edema and improves venous blood flow^{4,5}
- Inelastic compression can increase arterial perfusion in patients with mild to moderate peripheral arterial disease^{4,5}

Resting Pressure

Donning Compriflex



Resting Pressure

Laying Down Sitting Still

Working Pressure

Walking, Stair Climbing, Sit to Stand, Toe Raises

Safe, Precise Resting Pressure

- Inelastic compression can be used safely for patients with mild to moderate peripheral arterial disease^{4,5*}
- Accutab delivers precise dosages for controlled initial resting pressure when donning Sigvaris Compriflex garment



For more details,
please visit
sigvaris-accutab.com

Wrap Details

Variable Compression

Each Compreflex Calf wrap includes the Accutab, allowing users to accurately adjust the garment to the desired level of compression

Intelligent Design

Garment features a contoured calf, providing a superior fit

Premium Materials

Soft, conforming Breathe-O-Prene® provides a comfortable fit

Easy Application

With a front stretch panel to secure the garment in place and straps that roll back, the garment is easy to don and adjust, improving adherence

Superior Quality

Each garment is double stitched to ensure durability and is backed by a six-month warranty*



Basic Liner



Features

Non-compressive

- Easy to apply knee-high liners
- *drirelease*® moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties

Included With

- Compreflex Standard Calf & Foot
- Compreflex Standard Calf

Transition Liner



Features

10–15mmHg in Foot and Ankle Only

- Soft, flat seam knee-high liners
- *drirelease* moisture wicking and odor control
- Mild compression in foot and ankle only
- Padded sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

Included With

- Compreflex Transition

Complete Liner



Features

10–15mmHg in Foot, Ankle, and Calf

- Comfortable knee-high liners
- Mild compression in foot, ankle and calf
- *drirelease* moisture wicking and odor control
- Padded sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

Included With

- Compreflex Complete

Compreboot Lite



Features

- Ideal for mild edema/lymphedema
- Highly adjustable design (may be trimmed to fit)
- Ultra-thin profile with adjustable velcro tabs
- Breathable stretch fabric conforms to ankle
- Provides mild, therapeutic compression

Compreboot Standard



Features

- Ideal for mild-moderate edema/lymphedema
- Features neoprene material with velcro tabs
- Four-way stretch for firm compression around ankle
- Open heel and toe
- Provides moderate, therapeutic compression
- Included with Compreflex Standard Calf & Foot

Compreboot Plus



Features

- Ideal for chronic/acute edema/lymphedema
- Chipped foam lining
- Removable gripper patches provide traction
- Straps are removable for custom positioning or replacement
- Provides moderate, therapeutic compression

Strap Extender

Easy To Use

- Provides up to 10cm of additional circumferential size range to Compreflex calf wraps
- Potential inexpensive solution for patients outside of standard size ranges of Compreflex leg wraps, or needing custom size



Additional Compreflex Inelastic Wraps



Compreflex Standard Arm

Provides moderate, easy-to-adjust compression. Features inner sleeve for simple, one-handed application.



Compreflex Standard Knee

Soft, conforming compression.
Features stretch panel and
straps that close over patella.

*For more information on
our Compreflex line, visit
our website or call Customer
Care to request a 2019
SIGVARIS GROUP Catalog*

*PH: 800-322-7744
W: sigvaris.com*



Compreflex Standard Thigh

20–50mmHg of comfortable
compression. Includes hip
attachment and belt.



Compreflex Reduce System

Reduction system for the foot,
calf, knee and thigh. Patients
with lobules or non-conforming
limbs may use the lobe strap to
improve garment fit.

| Sizing Charts

Compreflex Calf Wraps – All Products (A-C1)

| | Small | Medium | Large | X Large | XX Large |
|----|----------------------------------|---------|---------|---------|----------|
| C1 | 29–39cm | 34–44cm | 40–50cm | 46–56cm | 54–64cm |
| C | 29–39cm | 34–44cm | 40–50cm | 46–56cm | 54–64cm |
| B | 24–34cm | 29–39cm | 34–44cm | 39–49cm | 44–55cm |
| A | 19–26cm | 21–30cm | 26–36cm | 31–41cm | 36–46cm |
| G | Regular: 30–33.9cm Tall: 34–42cm | | | | |

Compreboot (I–K)

| | Small | | Medium / Large | | X Large / XX Large | |
|---|---------|---------|----------------|----------|--------------------|----------|
| | Regular | Long | Regular | Long | Regular | Long |
| I | 28–36cm | 28–36cm | 39cm max | 39cm max | 44cm max | 44cm max |
| J | 22–26cm | 22–26cm | 30cm max | 30cm max | 33cm max | 33cm max |
| K | 14–18cm | 19–23cm | 18–20cm | 23–25cm | 20–22cm | 25–27cm |

Transition & Complete Liner

| Size | Small | Medium | Large | X Large | XX Large |
|--------------|---------|---------|---------|---------|----------|
| Men's Shoe | <6.5 | 6.5–8.5 | 9–11 | 11.5–14 | 11.5–14 |
| Women's Shoe | <7.5 | 7.5–9.5 | 10–12 | 12.5–15 | 12.5–15 |
| C | 29–39cm | 34–44cm | 40–50cm | 46–56cm | 54–64cm |
| A | 19–26cm | 21–30cm | 26–36cm | 31–41cm | 36–46cm |

| How To Order

Phone
800.322.7744, Ext. 1

Email
us_orders@sigvaris.com

Dealer Portal
sigvaris-online.com

Sigvaris Compreflex Order Form

| | | | | | |
|--------|-----------|---------|------|--|--|
| P.O. # | Account # | Contact | Date | | |
|--------|-----------|---------|------|--|--|

| | | | | |
|-------|--------------|-----|--------|--------|
| Phone | Patient Name | Age | Height | Weight |
|-------|--------------|-----|--------|--------|

Company Name

Dropship Name & Address

Product Information

| | | | |
|--|------|--------|---|
| Leg | Size | Length | Color |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | | <input type="checkbox"/> Black <input type="checkbox"/> Beige |

| | | | |
|--|------|--------|---|
| Leg | Size | Length | Color |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | | <input type="checkbox"/> Black <input type="checkbox"/> Beige |

Instructions

The measurements indicated are required when ordering Sigvaris compression garments. Complete this chart with the required measurements. Use the chart on page 10 to determine the optimal product size.

Product

- Compreflex Transition Calf**
Includes Transition Liner
- Compreflex Complete Calf**
Includes Complete Liner
- Compreflex Standard Calf & Foot**
Includes Compreboot
- Compreflex Standard Calf**
No foot piece included


Items Sold Separately

- Compreboot Lite**
- Compreboot Standard**
- Compreboot Plus**
Size: _____ Length: _____
- Transition Liner**
- Complete Liner**
- Basic Liner**
Size: _____ Length: _____
- Strap Extender**
Quantity: _____

Circumference

Left Right

| | | |
|----|-------|-------|
| C1 | _____ | _____ |
| C | _____ | _____ |
| B | _____ | _____ |
| A | _____ | _____ |
| I | _____ | _____ |
| J | _____ | _____ |
| K | _____ | _____ |



Customer Care

800.322.7744

Email

us_orders@sigvaris.com

Find a retailer

sigvaris.com/find

Dealer Portal

sigvaris-online.com



* Please consult your doctor, or other qualified health care professional before using any product discussed within the literature. See product insert for complete information.

(1) Partsch, Hugo, et al. "Measurement of lower leg compression in vivo: recommendations for the performance of measurements of interface pressure and stiffness." *Dermatologic Surgery* 32.2 (2006): 224-233.

(2) <http://www.tagungsmanagement.org/comp/ndex.php/lcc-consensus>

(3) Mosti, Giovanni, and Susanna Rossari. "The importance of measuring sub bandage pressure and presentation of new measuring device." *Acta Vulnol* 6 (2008): 31-36.

(4) Mosti, Giovanni, Maria Letizia Iabichella, and Hugo Partsch. "Compression therapy in mixed ulcers increases venous output and arterial perfusion." *Journal of vascular surgery* 55.1 (2012): 122-128.

(5) Mosti, Giovanni. "Compression in mixed ulcers: venous side" *Phlebology* 29.1_suppl (2014): 13-17.

** Garment includes a 6-month warranty against defects in materials and workmanship

Sigvaris, SIGVARIS GROUP, Compreflex, Comprefit, Compreboot, Cooflex, and Accutab are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide.

© Copyright 2019 SIGVARIS, Inc. 2019 drirelease. drirelease is an Optimizer Brands company. Breathe-O-Prene is a registered trademark of Accumed Innovative Technologies, LLC. All Rights Reserved.

LIT-COMPREFLEXBOOK 06.19