



People.Health.Care.

# TributeWrap™

Promoting skin health with foam technology.



**Easy-to-use finger grip**  
Simple pull function for the right, comfortable fit.

**Chevron channeling**  
Designed based on MLD principles of directing lymphatic fluid

**Durable VELCRO® brand fasteners**  
Robust hold for long-lasting wear.



People.Health.Care.

## TributeWrap™ Adjustable foam compression wraps.

TributeWrap™ is an adjustable garment with foam technology that nestles into the tissue to promote skin health. It comfortably adjusts to your patient's unique shape and lifestyle to help maintain gains made during therapy, while allowing the flexibility to support further reduction. Simply slide it on, adjust the straps, and cover with Sleep Sleeve for a comfortable experience.

-  Off-the-Shelf
-  Evening & Nighttime
-  Low Activity

### TributeWrap Glove

left or right orientation, sold individually

Size	Circumference	
	B <sup>°</sup>	C <sup>°</sup>
Small	16–20	16.5–21
Medium	18–22	17.5–22
Large	20–24.5	18.5–23

### TributeWrap Wrist to Axilla

left or right orientation, sold individually

Size	Length CG <sup>L</sup>	Circumference		
		C <sup>°</sup>	E <sup>°</sup>	G <sup>°</sup>
Small	43–47	16–21	22–28	25–31
Medium		17–22	26–32	30–37
Large		18–23	30–36	35–44

### TributeWrap Below Knee


left or right orientation, sold individually

Size	Length		Circumference	
	AY <sup>L</sup>	YD <sup>L</sup>	B <sup>°</sup>	C <sup>°</sup>
Small	19–22	40–44	24–30	30–36
Medium			27–33	33–41.5
Large			30–36	36–45

All measurements in centimeters.

### Indications

- Lymphedema
- Edema



**Proven Foam Technology**  
Breaks up and softens hardened or fibrotic tissue to support skin health.



**Sleep Sleeve NEW**  
Complimentary cover included with all styles to ensure the straps don't get tangled in the sheets.



POWERED BY   
Innovative interior fabric with chemical-free Cooling Power that rapidly wicks away sweat and moisture to regulate body temperature.



L&R USA INC. · 3880 W Wheelhouse Road, Milwaukee, WI 53208 USA · Phone: 855-892-4140 · Fax: 414-892-4150 · Email: inquiries@us.LRmed.com

Coolcore® is a registered trademark of Coolcore, LLC. VELCRO® is a registered trademark of Velcro BVBA.