

JOBST[®] Elvarex[®]

Custom-fit, flat-knit compression garments

Why JOBST®?

JOBST is the **#1 physician recommended brand** of medical compression legwear for the prevention, treatment, and management of venous and lymphatic conditions. The brand was founded by Conrad Jobst, mechanical engineer and industry pioneer, in 1950.



Innovation is at the heart of our product excellence. Product development efforts are directed towards delivering the highest quality product with the highest medical efficacy.

Our **product experts**, called Medical Education Advisors, and our devoted team of Sales Professionals, provide customized education and solutions for both fitters and therapists.

Knowledgeable customer and product service representatives provide a personalized touch to problem solving.

Quick and efficient processing from order entry to garment delivery. Elvarex® garments are delivered in an average of 5 business days. Two free of charge replacements on the initial order will be allotted within 45 days of confirmed delivery of parcel carrier for reasons other than manufacturing or material flaws.

The JOBST® Elvarex® Family



Each patient has different requirements and preferences which can be accommodated quickly and easily with the JOBST® Elvarex® custom-fit range. Custom-fit garments are manufactured to each patient's individual measurements so even patients with limb shape distortion are better treated as garments do not curl, twist, or cut into the skin.

Benefits of Flat-Knit Garments

The number of stitches can be varied to achieve virtually any desired shape to accommodate almost any level of limb distortion, providing:

- · Exact anatomical fit for each individual patient
- High breathability to enhance patient comfort
- Provides optimal compression profile, even when the edema is extreme
- Fabric does not dig in at skin folds
- High ratio of working pressure to resting pressure

JOBST[®] Elvarex[®]

The gold standard in custom, flat-knit compression garments

- Flat-knit structure bridges skin folds instead of working into them
- Latex core provides superior containment and allows garment to return to original shape for extended wear
- Excellent breathability for comfort and patient compliance
- Upper and lower limb garments for the treatment of moderate to severe lymphedema and chronic edema
- Provides SPF UV 80 protection



Compression Classes

Upper Extremities

CCL 1 = 15-21 mmHa* CCL 2 = 23-32 mmHg* CCL 2 Forte = 23-32 mmHg*

Lower Extremities

CCL 1 = 18-21 mmHg* CCL 2 = 23-32 mmHg* CCL 3 = 34-46 mmHg* CCL 3 Forte = 34-46 mmHg* $CCL 4 = 49-70 \text{ mmHg}^*$ CCL 4 Super = 60-90 mmHg*

Colors

Upper Extremities



Beige

Caramel

Cherry







Lower Extremities















NEW

Henna**1



Denim**1



Colored Seams

Liven up any lower extremity garment with a 'Colored Seam', a new option to add a bit of flare, style and individuality to our existing colors range.



Basic Styles



Sleeve



Gauntlet





Glove





Thiah



Pantyhose



Two-Lea





Chap style





Bermuda





Toe Cap

* Design Pressure ** Colored seams are not available in Henna, Denim and Graphite. † Not available in Toe Cap.

Knee

^{*} CAUTION: These products contains natural rubber latex which may cause allergic reactions.

JOBST® Elvarex® Soft

Gentle and latex-free

- Upper and lower limb garments for the treatment of mild to moderate lymphatic and venous conditions
- Designed for those patients with dry, sensitive or fragile skin
- Provides a more comfortable choice for patients
- Soft textured yarns create an extremely smooth fabric surface to aid ease of donning
- Combining skin-protecting softness and excellent breathability



Compression Classes

Upper Extremities

CCL 1 = 15-21 mmHg*

CCL 2 = 23-32 mmHg*

Lower Extremities

CCL 1 = 18-21 mmHg*

CCL 2 = 23-32 mmHg*

CCL 3 = 34-46 mmHg*

Colors

Upper Extremities



Beige

















Lower Extremities



















Colored Seams

Option to add a bit of flare, style and individuality to our lower extremity beige garment.



Basic Styles







Thigh











Bermuda

Chaps style

^{*} Design Pressure

JOBST® Elvarex® Plus

The cutting edge of flat-knit technology

- Seamless comfort that provides maximum range of motion without chafing or pressure points
- 3D knitting technology ensures an ideal anatomical fit that wears like a second skin
- Effective edema control with reliable therapeutic outcomes

Compression Classes

Upper Extremities

CCL 2 = 23-32 mmHg* CCL 1 = 15-21 mmHg*

Lower Extremities

CCL 1 = 18-21 mmHg* CCL 2 = 23-32 mmHg*



Colors

Upper Extremities



Beige



Black













Lower Extremities











Cherry





Basic Styles



AA Foot Cap



AC1 Glove with thumb and fingers to wrist



AC1 Gauntlet with thumb to wrist



AE Glove or gauntlet with thumb to elbow

JOBST® Elvarex® Soft Seamless

Soft and seamless, where it matters most

- · Soft, durable, and breathable flat-knitted fabric
- Flat edge finish at fingers and toes decreases bulk and prevents curling
- · Easy to don, comfortable to wear
- · Latex free

Compression Classes

Upper Extremities

CCL 1 = 15-21 mmHg* CCL 2 = 23-32 mmHg*

Lower Extremities

CCL 1 = 18-21 mmHg* CCL 2 = 23-32 mmHg*



Colors

Upper Extremities



Beige







Cherry











Lower Extremities

















Basic Styles



AA Foot Cap



AC1 Glove with thumb and fingers to wrist



AC1 Gauntlet with thumb to wrist



AE Glove or gauntlet with thumb to elbow

^{*} Design Pressure

JOBST® Elvarex® and JOBST® Elvarex® Soft Options

Functional Zone for greater wearing comfort

Give your patients more freedom! Additional Functional Zones can be knitted into the garments to give greater **mobility** and **support** around joints, according to the patient's needs.

The patented Functional Zone is a uniquely knitted two-layer construction created from selected soft and particularly elastic yarns. These reduce cutting in and wrinkling to provide the necessary **comfort** as well as retaining **compression**.



Functional zones available for:





Knee Comfort



Elbow Comfort



T-Heel



Pressure Panty



Slant Foot



Top Comfort



Top Bands

Stage 0 Stage 2 Stage 3 Stage 1 Latent: A subclinical state Mild: This represents early Moderate: Limb elevation Severe: The tissue is hard where swelling is not onset of the condition alone rarely reduces (fibrotic) and pitting is evident despite impaired where there is a accumulaabsent. Skin changes such swelling and pitting is lymph transport. This stage tion of tissue fluid that manifest. There may or as thickening, hypermay exist for months or subsides with limb may not be pitting as tissue pigmentation, increased skin years before edema elevation. This edema may fibrosis is more evident. folds, fat deposits and warty becomes evident. be pitting at this stage. overgrowths develop. Welt with SoftFit technology Knitted welt Silicone band



Seeking more information about vein diseases, lymphedema, and our products?

The JOBST® Compression Institute serves as a resource for healthcare providers and clinicians seeking further education as well as product information and resources. The JOBST® Compression Institute houses product literature, studies, order forms, as well as donning and meausring videos. Additionally, the JOBST® Compression Institute offers a variety of nationwide workshops and continuing education classes to train you on JOBST® products and proper measuring and fitting techniques. Up to 9 hours of CEU credits available. Register at www.jobstcompressioninstitute.com.

 $Please\ refer\ to\ the\ product\ label\ and\ /\ or\ package\ insert\ for\ full\ instructions\ on\ the\ safe\ use\ of\ these\ products.$





