

SIGVARIS  
GROUP

# 2019 Compression Solutions

Graduated Compression  
Socks & Hosiery  
for Venous-Lymphatic Disorders



### How to Order Medical Products

Series Number 23 = Cotton (See individual Catalog pages for series number)	Compression Level 1 = 15-20mmHg 2 = 20-30mmHg* 3 = 30-40mmHg 4 = 40-50mmHg** 5 = 50-60mmHg**	Product Style C = Knee-high or Calf (See individual Catalog pages for available styles)	Size ML = Medium Long (See individual Catalog pages for sizing)	Gender M = Men W = Women O = Open toe	Color 99 = Black (See individual Catalog pages for available colors)
<b>23</b>	<b>2</b>	<b>C</b>	<b>ML</b>	<b>M</b>	<b>99</b>

\* 18-23mmHg available in Anti-Embolism stockings only, 18-25mmHg available in Diabetic Compression Socks only

\*\* 40-50mmHg available in Secure stockings, 40-50mmHg and 50-60mmHg available in Natural Rubber only

**To order Plus Calf with Grip-top or Plus Size Pantyhose:**  
add "/P" after the color in the ordering code (860P for women only)

**To order Calf Style with Grip-top:**  
add "/S" after the color in the ordering code (see product Catalog pages for availability)

**To order Thigh with Waist Attachment:**  
add "/L" to indicate the Left leg or "/R" to indicate the Right leg after the color in the ordering code (860W and 500W only)

**To order Motion High Tech and Motion Thermoregulating Wool:**  
The letter for Gender is not used in the ordering code.

**Example:** 412CSM59 = Motion High Tech, 20-30mmHg, Closed toe Calf, Size Small Medium, Color Red

### How to Order WELL BEING & CONSUMER HEALTH Products

Series Number 120 = SHEER FASHION (See individual Catalog pages for series number)	Product Style N = Thigh-high (See individual Catalog pages for available styles)	Size A = Small (See individual Catalog pages for sizing)	Color 33 = Natural (See individual Catalog pages for available colors)
<b>120</b>	<b>N</b>	<b>A</b>	<b>33</b>

**To order ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, and MERINO OUTDOOR SOCKS:**

Size is only one letter and the letter for Gender is not used in the ordering code.

**Example:** 421CL31 = MERINO OUTDOOR SOCKS, 15-20mmHg, Closed toe Calf, Size Large, Color Olive

**To Order:**  
1.800.322.7744, ext. 1

**Online Dealer Portal:**  
[sigvaris-online.com](http://sigvaris-online.com)

Previous	Series Number	New	Segment	
COTTON	230	Cotton	Essential	
SELECT COMFORT	860	Opaque		
SEA ISLAND COTTON	220	Sea Island Cotton	Style	
ALL-SEASON MERINO WOOL	240	Merino Wool		
ALLURE	710	Patterns		
MIDSHEER	750	Medium Sheer		
EVERSHEER	780	Sheer		
MIDTOWN MICROFIBER	820	Microfiber		
MICROFIBER SHADES	830	Patterns		
SOFT OPAQUE	840	Soft Opaque		
CUSHIONED COTTON	360	Cushioned Cotton		Motion
TRAVERSE SOCKS	412	High Tech		
MERINO WOOL	422	Thermoregulating Wool		
DAILY COMFORT	850	Comfort		

Previous	New	Category
COMPREFLEX LITE TRANSITION	Compreflex Transition Calf	Medically Complex Edema
COMPREFLEX LITE COMPLETE	Compreflex Complete Calf	
COMPREFLEX BK	Compreflex Standard Calf & Foot	
COMPREFLEX NF	Compreflex Standard Calf	
COOLFLEX NF	Coolflex Standard Calf	

**Please follow these instructions carefully to determine the correct size for different product styles.**

Measurements should always be taken in the following order, and as soon after waking as possible when swelling is minimal.

### CALF STYLE

Open/Closed Toe

<b>1</b>	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. <b>Note: This should be your first measurement.</b>
<b>2</b>	Measure the calf circumference at the fullest part of the calf.
<b>3</b>	Determine the calf length by measuring from the floor up to the fibular head. This measurement is taken on the outside of the leg.

### PANTYHOSE AND THIGH-HIGH STYLE

Open/Closed Toe


<b>1</b>	Measure the ankle circumference directly above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE, or X-LARGE. <b>Note: This should be your first measurement.</b>
<b>2</b>	Measure the calf circumference at the fullest part of the calf.
<b>3</b>	Skip Step 3 (but refer to General Guidelines #6 below).
<b>4</b>	Measure the thigh at its widest circumference.
<b>5</b>	Determine the leg length by measuring from the floor up to the gluteal fold.
<b>6</b>	Determine the hip circumference by measuring the widest part of the hips.

### GENERAL GUIDELINES

<b>1</b>	Measurements should always be taken as soon after waking as possible, when swelling is minimal.
<b>2</b>	Measure against bare skin, not over clothing.
<b>3</b>	The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to the right for examples).
<b>4</b>	Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
<b>5</b>	Length measurements should be taken from the ground up.
<b>6</b>	Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles.




**1.** Ankle circumference directly above ankle bone



**2.** Calf circumference at fullest part of the calf



**3.** Calf length from the floor to the fibular head



**4.** Widest circumference of the thigh



**5.** Leg length from the floor up to the gluteal fold



**6.** Circumference at the hip





SIGVARIS  
GROUP

Established in 1864, SIGVARIS GROUP is, and has always been, a family-owned company. Founded and headquartered in Switzerland, our passion for helping people feel their best is why we've been a global market leader for 155 years and counting.

Our products restore, assist and enhance people's lives through brands like:

**sigvaris**

Medical Compression Socks, Hosiery  
and Solutions for Medically Complex Edema

**SIGVARIS WELL BEING**

Consumer Healthcare: Simplified Sizing  
Compression Socks and Hosiery

We are renowned for making exceptional products that combine advanced technology and quality craftsmanship with diverse styles that fit our modern customers' lives.

When heritage, competence and dedication come together, there's no limit to where we'll go.

**Live your life. Feel your best.**

For more information, please visit:  
**[sigvaris.com](https://sigvaris.com)**

# sigvaris

	SERIES	PAGE	COMPRESSION LEVELS	STYLES
Essential	Cotton <b>(230)</b>	8	20–30mmHg 30–40mmHg	<b>closed toe:</b> calf, calf w/grip-top, thigh-high w/grip-top <b>open toe:</b> calf, thigh-high w/grip-top
	Opaque <b>(860)</b>	10	20–30mmHg 30–40mmHg	<b>closed toe:</b> calf, petite calf, calf plus, calf w/grip-top, thigh w/grip-top, petite thigh-high w/grip-top, extra-large thigh w/grip-top, pantyhose, petite pantyhose, plus-sized pantyhose <b>open toe:</b> calf, calf plus, calf w/grip-top, thigh w/grip-top, extra-large thigh w/grip-top, thigh-high w/waist attachment
Style	Sea Island Cotton <b>(220)</b>	ZÜRICH COLLECTION 14	20–30mmHg	<b>closed toe:</b> calf
	Merino Wool <b>(240)</b>	ZÜRICH COLLECTION 16	20–30mmHg	<b>closed toe:</b> calf
	Patterns <b>(710)</b> <i>(for women only)</i>	18	15–20mmHg 20–30mmHg	<b>closed toe:</b> thigh-high w/grip-top, pantyhose
	Medium Sheer <b>(750)</b> <i>(for women only)</i>	20	20–30mmHg	<b>closed toe:</b> calf, thigh-high w/grip-top, pantyhose
	Sheer <b>(780)</b> <i>(for women only)</i>	22	15–20mmHg 20–30mmHg 30–40mmHg	<b>closed &amp; open toe:</b> calf, thigh-high w/grip-top, pantyhose, maternity pantyhose
	Microfiber <b>(820)</b> <i>(for men only)</i>	24	15–20mmHg 20–30mmHg 30–40mmHg	<b>closed toe:</b> calf, calf w/grip-top, thigh w/grip-top, extended foot sizes
	Patterns <b>(830)</b>	26	20–30mmHg	<b>closed toe:</b> calf
	Soft Opaque <b>(840)</b> <i>(for women only)</i>	28	15–20mmHg 20–30mmHg 30–40mmHg	<b>closed toe:</b> calf, calf w/grip-top, thigh-high w/grip-top, pantyhose, maternity pantyhose <b>open toe:</b> calf, thigh-high w/grip-top, pantyhose
Motion	Cushioned Cotton <b>(360)</b>	32	20–30mmHg	<b>closed toe:</b> calf
	High Tech <b>(412)</b>	34	20–30mmHg	<b>closed toe:</b> calf
	Thermoregulating Wool <b>(422)</b>	36	20–30mmHg	<b>closed toe:</b> calf
	Comfort <b>(850)</b>	38	20–30mmHg	<b>closed toe:</b> calf
Specialty	Natural Rubber <b>(500)</b>	42	30–40mmHg 40–50mmHg 50–60mmHg	<b>open toe:</b> calf, thigh-high, thigh w/grip-top, thigh-high w/waist attachment, pantyhose
	Secure <b>(550)</b>	44	20–30mmHg 30–40mmHg 40–50mmHg	<b>closed &amp; open toe:</b> calf w/grip-top, thigh-high w/ grip-top
	Secure Armsleeve, Glove, and Gauntlet <b>(560)</b>	46	15–20mmHg 20–30mmHg 30–40mmHg	armsleeve w/ grip-top (glove and gauntlet available separately)
	Diabetic Compression Socks <b>(602)</b>	48	18–25mmHg	<b>closed toe:</b> calf
	Advance Armsleeve <b>(910)</b>	50	15–20mmHg 20–30mmHg 30–40mmHg	armsleeve w/ grip-top w/ gauntlet, armsleeve w/ grip-top w/out gauntlet
	Anti-Embolism Stockings <b>(930)</b>	52	18–23mmHg	<b>open toe:</b> calf, thigh-high

# SIGVARIS WELL BEING & CONSUMER HEALTH

SERIES	PAGE	COMPRESSION LEVEL	STYLES
CASUAL COTTON (146/186)	56	15–20mmHg	closed toe: calf
CUSHIONED COTTON (142/182)	58	15–20mmHg	closed toe: calf
SEA ISLAND COTTON (151/191) <small>ZÜRICH COLLECTION</small>	60	15–20mmHg	closed toe: calf
ALL-SEASON MERINO WOOL (152/192) <small>ZÜRICH COLLECTION</small>	62	15–20mmHg	closed toe: calf
MICROFIBER SHADES (143/183)	64	15–20mmHg	closed toe: calf
SHEER FASHION (120) <i>(for women only)</i>	66	15–20mmHg	closed toe: calf, thigh-high w/grip-top, pantyhose, maternity pantyhose open toe: calf
SOFT SILHOUETTE LEGGINGS (170) <i>(for women only)</i>	68	15–20mmHg	open foot: leggings
BUSINESS CASUAL (189) <i>(for men only)</i>	70	15–20mmHg	closed toe: calf
ATHLETIC RECOVERY SOCKS (401)	72	15–20mmHg	closed toe: calf
PERFORMANCE SLEEVES (412)	74	20–30mmHg	open foot: calf sleeve
MERINO OUTDOOR SOCKS (421)	76	15–20mmHg	closed toe: calf
EVERSOFT DIABETIC SOCKS (160)	78	8–15mmHg	closed toe: calf

Medically Complex Edema	82
Accessories	92

Custom sizes are available for a cotton product, a microfiber product, and our Natural Rubber series. For more information, please speak with your SIGVARIS GROUP Sales Representative or Customer Care:

**1.800.322.7744, ext. 1.**

# sigvaris

---

## Life changes. Your lifestyle doesn't have to.

At Sigvaris, we care about people. For 155 years, our compression products have offered effective and reliable solutions for venous conditions – anywhere, anytime, and for anyone.

As a Sigvaris Certified Fitter, when you provide **the Sigvaris solution** to your clients, you're helping them move with life, without compromise. Because you care about people too. And that's why we're proud to work with Fitters like you.

Our shared passion for helping people is what got us here, and we're excited to take our next step towards wellness together.

Let's take a look at what we've got in store.



# Be yourself.

## ESSENTIAL

### Support Yourself

Sigvaris Essential solutions always meet your customers' needs. These straightforward stockings are an everyday go-to for people who seek the comfort, support and reliability of traditional therapeutic socks and stockings.

#### Why choose **Essential**?

- Reduce swelling and discomfort.
- Long-distance travel.
- Everyday comfort and reliability.

## STYLE

### Express Yourself

Choice is empowering, and fashion can be as unique as you are. With Sigvaris Style, personality will shine as your customers select from distinctive hosiery and socks.

#### Why choose **Style**?

- Fashionable, modern patterns and designs.
- Special occasion and formal wear.
- Our most luxurious feel.

## MOTION

### Push Yourself

When daily comfort for hard-working legs is required, put your customers in motion. Sigvaris Motion socks offer **all day strong** support for people with active lives.

#### Why choose **Motion**?

- Cushioned sole for added comfort.
- Durable thermoregulating wool.
- Moisture-wicking yarns.



# sigvaris





## ESSENTIAL

### **Support** Yourself

Why choose **Essential**?

- Reduce swelling and discomfort
- Long-distance travel
- Everyday comfort and reliability

#### Essential Products

Previous	Series Number	New	Page
COTTON	<b>230</b>	Cotton	8
SELECT COMFORT	<b>860</b>	Opaque	10

230

## Essential COTTON

for women & men

Essential Cotton combines all-day cotton comfort and a classic ribbed design that is ideal for sensitive skin.

### Compression Level

20-30mmHg

30-40mmHg

### Features and Benefits



**Sensitive  
Skin**



**Cotton**



**Breathable**

- Supima® cotton for highest level of comfort and breathability
- Ideal for sensitive skin and allergies
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

### Fiber Content


- 50% Nylon, 25% Spandex, 25% Supima Cotton (content can vary slightly by style)

- Non-latex

**SUPIMA®**  
WORLD'S FINEST COTTONS

### Best Suited for

- Containment of various types of edema, including venous edema, phlebolympheidema, and early stages of lymphedema
- Ideal for dress or casual shoes
- Ideal for long-distance travel



**Silicone**  
band grip-top  
in calf and thigh

20-30mmHg  
Closed Toe Calf w/ Grip-top  
Medium Long in Black

**232CMLM99/S**



20-30mmHg  
Closed Toe Calf Medium Short in Black Mist



**232CMSW14**

**Sizing Chart**

	Small		Medium		Large		X Large	
<b>Ankle Circumference</b>	7"-8.5" (18-21.5cm)		8.5"-10" (21.5-25.5cm)		10"-11.5" (25.5-29cm)		11.5"-14" (29-36cm)	
<b>Calf Circumference</b>	11"-15.5" (28-39cm)		12.5"-17.5" (31.5-45cm)		14"-20" (35.5-51cm)		16"-24" (41-61cm)	
<b>Calf Length</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
<b>Thigh Circumference</b>	18"-26" (45.5-66cm)		19"-28.5" (48-73cm)		20.5"-33" (52cm-84cm)		NA	
<b>Leg Length: Women &amp; Open Toe</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>		
	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)		
<b>Leg Length: Men</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>		
	≤30" (76cm)	>30" (76cm)	≤30" (76cm)	>30" (76cm)	≤30" (76cm)	>30" (76cm)		
<b>Size</b>	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12 Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

**Styles and Colors**

<b>Closed Toe = ●</b> <b>Open Toe = ○</b>						
	WHITE 00	NAVY 10	BLACK MIST 14	CRISPA 66	CHOCOLATE 88	BLACK 99
<b>Calf (C)</b>	○	●	●	●○	●*	●
<b>Calf w/Grip-Top (C/S)</b> Silicone band grip-top				●*		●**
<b>Thigh w/Grip-Top (N)</b> Silicone band grip-top				●○		●

\* 20-30mmHg only

\*\* 30-40mmHg calf with grip-top only available in black, 30-40mmHg extra large calf only available in black with grip-top

Product Code Example: **232NMLW66**

Series	mmHg Level	Style	Size	Gender/ Open Toe	Color
23	2	N	ML	W	66

Product Code Example: **232CLSM66/S**

Series	mmHg Level	Style	Size	Gender/ Open Toe	Color	Grip-top
23	2	C	LS	M	66	/S

## Essential

### OPAQUE

for women & men

Essential Opaque offers durability you can trust in our largest range of sizes, making it your go-to choice for everyday hosiery.

### Compression Level


20-30mmHg

30-40mmHg

### Features and Benefits



**Durable**



**Largest  
Size Range**



**Opaque**


- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Cotton-lined gusset for pantyhose styles
- Plus size pantyhose great alternative for maternity

### Fiber Content

- 64% Nylon, 36% Spandex  
(content can vary slightly by style)
- Non-latex

### Best Suited for

- Containment of various types of edema including venous edema, phlebolympheidema, and early stages of lymphedema
- Post-vein treatment including thermal and non-thermal techniques
- Patients looking for petite and plus-sized options



**Silicone band**  
band grip-top  
in thigh and  
calf w/ Grip-top

20-30mmHg  
Open Toe Thigh Medium Long in Crispa

**862NML066**



30-40mmHg  
Closed Toe Calf Small Long in Black


**863CSLW99**

**Sizing Chart**

	Small			Medium			Large		X Large	
<b>Ankle Circumference</b>	7"-8.5" (18-21.5cm)			8.5"-10" (21.5-25.5cm)			10"-11.5" (25.5-29cm)		11.5"-14" (29-36cm)	
<b>Calf Circumference</b>	11"-15.5" (28-39cm)			12.5"-17.5" (31.5-45cm)			14"-20" (35.5-51cm)		16"-24" (41-61cm)	
<b>Plus Size Calf Circumference</b>	NA			15.5"-21" (39.5-53cm)			18"-24" (45.5-61cm)		20"-28" (51cm-71cm)	
<b>Calf Length</b>	<b>Petite*</b>	<b>Short</b>	<b>Long</b>	<b>Petite*</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>
	≤14" (35.5cm)	≤16" (41cm)	>16" (41cm)	≤16" (35.5cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
<b>Thigh Circumference</b>	18"-26" (45.5-66cm)			19"-28.5" (48-73cm)			20.5"-33" (52cm-84cm)		25"-37.5" (64cm-95cm)	
<b>Leg Length: Women &amp; open Toe</b>	<b>Petite*</b>	<b>Short</b>	<b>Long</b>	<b>Petite*</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>
	≤25.5" (65cm)	≤28.5" (72cm)	>28.5" (72cm)	≤25.5" (65cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)
<b>Leg Length: Men</b>	<b>Short</b>		<b>Long</b>	<b>Short</b>		<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>
	≤30" (76cm)		>30" (76cm)	≤30" (76cm)		>30" (76cm)	≤30" (76cm)	>30" (76cm)	≤30" (76cm)	>30" (76cm)
<b>Hip Circumference Guidelines</b>	32.5"-44.5" (83-112.5cm)			34.5"-50" (87-127cm)			38.5"-59.5" (98-150.5cm)		NA	
<b>Size</b>	SP*	SS	SL	MP*	MS	ML	LS	LL	XS	XL
<b>Plus Size Calf</b>					MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

**Styles and Colors**

Closed Toe = ●  
Open Toe = ○



<b>Women's Calf (C)</b>	●	●	○	●○
<b>Men's Calf (C)</b>			●	●
<b>Calf Plus w/Grip-Top (C/P)</b> <small>Only available with silicone band grip-top</small>			●○	●
<b>Calf w/Grip-Top (C/S)</b> <small>Silicone band grip-top</small>			●○	●
<b>Women's &amp; Open Toe Thigh w/Grip-Top (N)</b> <small>Silicone band grip-top</small>	●	●	○	●
<b>Men's Thigh w/Grip-Top (N)</b> <small>Silicone band grip-top</small>			●	●
<b>Pantyhose (P)</b>	●	●		●
<b>Pantyhose Plus (P/P)</b>			●	●
<b>Thigh w/Waist Attachment (W)</b> <small>Indicate leg LEFT (/L) or RIGHT (/R)</small>			○	

**Women's Shoe Sizes:** Petite Leg 4-8, Short Leg 5.5-9.5, Long Leg 8-12

**Men's Shoe Sizes:** Short Leg 7-11, Long Leg 9.5-14

**Plus (Maternity) Hip Circumference Guidelines** Small ≤56" (142cm), Medium ≤62" (157.5cm), Large ≤70" (178cm)

**\* Petite sizing (SP & MP):** available in women's closed toe only; without grip-top, 20-30mmHg only; Natural, Suntan, and Black only; unavailable in Plus sizes

**Product Code Example: 863PMSW33**

Series	mmHg Level	Style	Size	Gender/Open Toe	Color
86	3	P	MS	W	33

**Product Code Example: 862CLLW66/P**

Series	mmHg Level	Style	Size	Gender/Open Toe	Color	Grip-top
86	2	C	LL	W	66	/P



# sigvaris





## STYLE

### Express Yourself

Why choose **Style**?

- Fashionable, modern patterns and designs
- Special occasion and formal wear
- Our most luxurious feel

#### Style Products

Previous	Series Number	New	Page
SEA ISLAND COTTON	<b>220</b>	Sea Island Cotton	14
ALL-SEASON MERINO WOOL	<b>240</b>	Merino Wool	16
ALLURE	<b>710</b>	Patterns	18
MIDSHEER	<b>750</b>	Medium Sheer	20
EVERSHEER	<b>780</b>	Sheer	22
MIDTOWN MICROFIBER	<b>820</b>	Microfiber	24
MICROFIBER SHADES	<b>830</b>	Patterns	26
SOFT OPAQUE	<b>840</b>	Soft Opaque	28

220

## Style

### SEA ISLAND COTTON

for women & men

Style Sea Island Cotton socks are made with the world's rarest cotton, ensuring a luxurious feel that is as soft as cashmere.

### Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 60

### Features and Benefits



**Soft**



**Cotton**



**Breathable**

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

### Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex  
7% Supima Cotton
- Non-latex

### Best Suited for

- Management of symptoms related to chronic venous disease
- Patients who prefer natural fibers for everyday wear
- Long-distance travel



20-30mmHg  
Closed Toe Calf Small Long in Brown

**222CSLW11**



20-30mmHg  
Closed Toe Calf Medium Long in Black

**222CMLM99**




Sizing chart

	Small		Medium		Large		X Large	
Ankle Circumference	7"-8.5" (18-21.5cm)		8.5"-10" (21.5-25.5cm)		10"-11.5" (25.5-29cm)		11.5"-14" (29-36cm)	
Calf Circumference	11"-15.5" (28-39cm)		12.5"-17.5" (31.5-45cm)		14"-20" (35.5-51cm)		16"-24" (41-61cm)	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12    Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

Styles and colors

Closed Toe = ●

	 NAVY 10	 BROWN 11	 BLACK 99
Calf (C)	●	●	●

Product Code Example: **222CLLM99**

Series	mmHg Level	Style	Size	Gender	Color
22	2	C	LL	M	99

240

## Style

### MERINO WOOL

for women & men

Style Merino Wool socks combine extra-fine Australian Merino wool with a stylish design, making it the ultimate sock for year-round wear.

#### Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 62

#### Features and Benefits



Wool



Thermoregulating



Stylish



- Made with the highest quality extra-fine Australian Merino wool
- Thermoregulating for year-round wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

#### Fiber Content

- 53% Merino Wool, 38% Nylon, 9% Spandex
- Non-latex

#### Best Suited for

- Management of symptoms related to chronic venous disease
- Patients who prefer fashionable and premium fabric options
- Long-distance travel



20-30mmHg  
Closed Toe Calf Medium Long in Brown

**242CMLM11**



20-30mmHg  
Closed Toe Calf Medium Long in Navy

**242CMLW10**


**Sizing Chart**

	Small		Medium		Large		X Large	
Ankle Circumference	7"-8.5" (18-21.5cm)		8.5"-10" (21.5-25.5cm)		10"-11.5" (25.5-29cm)		11.5"-14" (29-36cm)	
Calf Circumference	11"-15.5" (28-39cm)		12.5"-17.5" (31.5-45cm)		14"-20" (35.5-51cm)		16"-24" (41-61cm)	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12    Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

**Styles and Colors**

Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---

Product Code Example: **242CMSW11**

Series	mmHg Level	Style	Size	Gender	Color
24	2	C	MS	W	11

710

## Style PATTERNS

for women

Style Pattern hosiery is available in fashionable patterns, making it a staple for your professional and everyday wear.

### Compression Level

15-20mmHg

20-30mmHg

### Features and Benefits



**Sheer**



**Elegant**



**Breathable**

- Stylish diamond-dot and stripe patterns for professional and dress wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Cotton-lined gusset for pantyhose styles

### Fiber Content

- 64% Nylon , 36% Spandex  
(Content can vary slightly by style)
- Non-latex

### Best Suited for


- Management of early stages of chronic venous disease
- Pre- and post-vein treatment with minimal edema
- Patients who prefer fashionable options
- Ideal for business or dress



**NEW!**




Scalloped Sensinnov® Thigh Band on Stripe patterns



**Sensinnov** thigh band on Diamond Dot patterns

15-20mmHg  
Closed Toe Thigh Medium Long in Graphite

**711NMLW91**



20-30mmHg  
Closed Toe Pantyhose Small Long in Onyx Stripe

**712PSLW19**








**Sizing Chart**

	Small		Medium		Large	
Ankle Circumference	7"-8.5" (18-21.5cm)		8.5"-10" (21.5-25.5cm)		10"-11.5" (25.5-29cm)	
Calf Circumference	11"-15.5" (28-39cm)		12.5"-17.5" (31.5-45cm)		14"-20" (35.5-51cm)	
Thigh Circumference	18"-26" (45.5-66cm)		19"-28.5" (48-73cm)		20.5"-33" (52cm-84cm)	
Leg Length	Short	Long	Short	Long	Short	Long
	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)
Hip Circumference Guidelines	32.5"-44.5" (83-112.5cm)		34.5"-50" (87-127cm)		38.5"-59.5" (98-150.5cm)	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes Short leg 5.5-9.5, Long leg 8-12

**Styles and colors**

Closed Toe = ●					
	DARK NAVY STRIPE 13	GRAPHITE STRIPE 17	ONYX STRIPE 19	GRAPHITE 91	BLACK 99
Thigh w/Grip-Top (N) Sensinnov band grip-top	●	●	●	●	●
Pantyhose (P)	●	●	●	●	●

Product Code Example: **712NLSW17**

Series	mmHg Level	Style	Size	Gender	Color
71	2	N	LS	W	17

750

## Style

### MEDIUM SHEER

for women

Thanks to the durability of Style Medium Sheer, you will have the confidence to wear sheer hosiery every day.

#### Compression Level

20-30mmHg

#### Features and Benefits



**Medium Sheer**



**Durable**



**Breathable**

- Fashionable, durable, medium sheer hosiery for everyday wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
  - Cotton-lined gusset for pantyhose styles

#### Fiber Content

- 68% Nylon, 32% Spandex (Content can vary slightly by style)
- Non-latex

#### Best Suited for

- Management of early stages of chronic venous disease
- Pre- and post-vein treatment with minimal edema
- Daily wear for business or casual
- Active professionals needing a durable sheer



**Silicone**  
band grip-top  
in thigh

20-30mmHg  
Closed Toe Thigh Medium Long in Black

**752NMLW99**



20-30mmHg  
Closed Toe Calf Small Long in Nightshade

**752CSLW94**

**Sizing chart**

	Small		Medium		Large	
Ankle Circumference	<b>7"-8.5"</b> (18-21.5cm)		<b>8.5"-10"</b> (21.5-25.5cm)		<b>10"-11.5"</b> (25.5-29cm)	
Calf Circumference	<b>11"-15.5"</b> (28-39cm)		<b>12.5"-17.5"</b> (31.5-45cm)		<b>14"-20"</b> (35.5-51cm)	
Calf Length	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
Thigh Circumference	<b>18"-26"</b> (45.5-66cm)		<b>19"-28.5"</b> (48-73cm)		<b>20.5"-33"</b> (52cm-84cm)	
Leg Length	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>
	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)
Hip Circumference Guidelines	<b>32.5"-44.5"</b> (83-112.5cm)		<b>34.5"-50"</b> (87-127cm)		<b>38.5"-59.5"</b> (98-150.5cm)	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes Short 5.5-9.5, Long 8-12

**Styles and colors**

Closed Toe = ●



	DARK NAVY 08	NATURAL 33	NIGHT-SHADE 94	BLACK 99
Calf (C)	●	●	●	●
Thigh w/Grip-Top (N) <small>Silicone band grip-top</small>	●	●	●	●
Pantyhose (P)	●	●	●	●

Product Code Example: **752NMLW33**

Series	mmHg Level	Style	Size	Gender	Color
75	2	N	ML	W	33

780

## Style

### **SHEER**

#### for women

Style Sheer hosiery offers fashionable sheerness you can feel comfortable in, making it ideal for dress and professional wear.

### Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

### Features and Benefits



**Sheer**



**Elegant**



**Breathable**


- Fashionable for professional and dress wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Cotton-lined gusset for pantyhose styles

### Fiber Content

- 15-20, 20-30mmHg: 69% Nylon , 31% Spandex (content can vary slightly by style)
- 30-40mmHg: 60% Nylon , 40% Spandex (content can vary slightly by style)
- Non-latex

### Best Suited for

- Management of early stages of chronic venous disease
- Pre- and post-vein treatment with minimal edema
- Patients who prefer fashionable options
- Ideal for business or dress



**Sensinnov**  
band grip-top  
in thigh

15-20mmHg  
Closed Toe Thigh Medium Long in Dark Navy

**781NMLW08**



20-30mmHg  
Closed Toe Calf Medium Long in Mocha

**782CMLW85**

**Sizing chart**

	Small		Medium		Large	
Ankle Circumference	7"-8.5" (18-21.5cm)		8.5"-10" (21.5-25.5cm)		10"-11.5" (25.5-29cm)	
Calf Circumference	11"-15.5" (28-39cm)		12.5"-17.5" (31.5-45cm)		14"-20" (35.5-51cm)	
Calf Length	Short	Long	Short	Long	Short	Long
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
Thigh Circumference	18"-26" (45.5-66cm)		19"-28.5" (48-73cm)		20.5"-33" (52cm-84cm)	
Leg Length	Short	Long	Short	Long	Short	Long
	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)
Hip Circumference Guidelines	32.5"-44.5" (83-112.5cm)		34.5"-50" (87-127cm)		38.5"-59.5" (98-150.5cm)	
Size	SS	SL	MS	ML	LS	LL

**Women's Shoe Sizes** Short 5.5-9.5, Long 8-12

**Maternity Hip Circumference Guidelines** Small ≤56" (142cm), Medium ≤62" (157.5cm), Large ≤70" (178cm)

**Styles and Colors**

	Closed Toe = ● Open Toe = ○								
Calf (C)	●	●○	●○	●○	●○	●○	●○	●	●○
Thigh w/Grip-Top (N) <small>Sensinnov band grip-top</small>	●	●○	●○	●○	●○	●○	●○	●	●○
Pantyhose (P)	●	●○	●○	●○	●○	●○	●○	●	●○
Maternity Pantyhose (M)**		●							●

\* 15-20mmHg & 20-30mmHg ONLY

\*\* 20-30mmHg ONLY

**Product Code Example: 781NMLW36**

Series	mmHg Level	Style	Size	Gender/Open Toe	Color
78	1	N	ML	W	36

**Product Code Example: 783PSSO33**

Series	mmHg Level	Style	Size	Gender/Open Toe	Color
78	3	P	SS	O	33

820

## Style

### MICROFIBER

for men

Style Microfiber socks offer microfiber softness for your professional and everyday wear.

### Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

### Features and Benefits



**Microfiber**



**Soft**



**Stylish**


- 100% of microfiber sits next to skin for a soft feel
- Excellent breathability and moisture management for added comfort
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

### Fiber Content

- 87% Nylon , 13% Spandex  
(content can vary slightly by style)
- Non-latex

### Best Suited for

- Management of early stages of chronic venous disease
- Patients who prefer fashionable options
- Men with larger feet
- Long-distance travel



**Silicone grip-**  
top band in calf\*

15-20mmHg  
Closed Toe Calf Medium Long in Navy

**821CMLM10**



**Sensinno**  
band grip-top  
in thigh

30-40mmHg  
Closed Toe Thigh Medium Long in Tan Khaki

**822NMLM32**



**Sizing Chart**

	Small		Medium		Large		X Large	
<b>Ankle Circumference</b>	<b>7"-8.5"</b> (18-21.5cm)		<b>8.5"-10"</b> (21.5-25.5cm)		<b>10"-11.5"</b> (25.5-29cm)		<b>11.5"-14"</b> (29-36cm)	
<b>Calf Circumference</b>	<b>11"-15.5"</b> (28-39cm)		<b>12.5"-17.5"</b> (31.5-45cm)		<b>14"-20"</b> (35.5-51cm)		<b>16"-24"</b> (41-61cm)	
<b>Calf Length</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
<b>Thigh Circumference</b>	<b>18"-26"</b> (45.5-66cm)		<b>19"-28.5"</b> (48-73cm)		<b>20.5"-33"</b> (52cm-84cm)		<b>NA</b>	
<b>Leg Length</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>		
	≤30" (76cm)	>30" (76cm)	≤30" (76cm)	>30" (76cm)	≤30" (76cm)	>30" (76cm)		
<b>Size Shoe size</b>	<b>SS</b> 7-11	<b>SL</b> 9.5-14	<b>MS</b> 7-11	<b>ML</b> 9.5-14	<b>LS</b> 7-11	<b>LL</b> 9.5-14	<b>XS</b> 7-11	<b>XL</b> 9.5-14
				<b>MX*</b> 12.5-17.5		<b>LX*</b> 12.5-17.5		<b>XX*</b> 12.5-17.5


**Men's shoe sizes:** Short leg 7-11, long leg 9.5-14 (For extended foot sizes: long leg 12.5-17.5)

\* MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15-20 & 20-30mmHg only.

Extended foot sizes not available with grip top

**Styles and Colors**

Closed Toe = ●



	NAVY 10***	TAN KHAKE 32	STEEL GREY 93***	BLACK 99
<b>Calf (C)</b>	●	●	●	●
<b>Calf w/Grip-Top (C/S)**</b> <small>Silicone band grip-top</small>		●		●
<b>Thigh w/Grip-Top (N)</b> <small>Sensinnov band grip-top</small>		●		●

\*\* Calf w/ Grip-Top not available in 15-20mmHg

\*\*\* Navy and Steel Grey colors only available in 15-20mmHg and 20-30mmHg

**Product Code Example: 822NLSM99**

Series	mmHg Level	Style	Size	Gender	Color
82	2	N	LS	M	99

**Product Code Example: 822CMLM32/S**

Series	mmHg Level	Style	Size	Gender	Color	Grip-top
82	2	C	ML	M	32	/S

830

## Style PATTERNS

for women & men

Style Microfiber Patterns socks combine ultimate comfort with a variety of fashionable colors and designs to spice up your wardrobe.

### Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 62

### Features and Benefits



**Microfiber**



**Stylish**



**Patterns**

- 100% of microfiber sits next to skin for a soft feel
- Fashionable, durable socks for everyday wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

### Fiber Content

- 87% Nylon, 13% Spandex  
(content can vary slightly by style)
- Non-latex

### Best Suited for

- Management of symptoms related to chronic venous disease
- Patients who prefer fashionable options
- Long-distance travel



20-30mmHg  
Closed Toe Calf Medium Short in Pink Stripe

**832CMSW15**



20-30mmHg  
Closed Toe Calf Large Long in Graphite Heather

**832CLLM27**

**Sizing Chart**

	Small		Medium		Large	
Ankle Circumference	7"-8.5" (18-21.5cm)		8.5"-10" (21.5-25.5cm)		10"-11.5" (25.5-29cm)	
Calf Circumference	11"-15.5" (28-39cm)		12.5"-17.5" (31.5-45cm)		14"-20" (35.5-51cm)	
Calf Length	Short	Long	Short	Long	Short	Long
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12    Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

**Styles and Colors**

Closed Toe = ●								
Calf (C)	●	●	●	●	●	●	●	●
Closed Toe = ●								
Calf (C)	●	●	●	●	●	●		

Product Code Example: **832CMLW48**

Series	mmHg Level	Style	Size	Gender	Color
83	2	C	ML	W	48

## Style **SOFT OPAQUE** for women

Style Soft Opaque hosiery blends a fashionable opaque look with long-lasting microfiber softness for your professional and everyday wear.

### Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

### Features and Benefits



**Opaque**



**Microfiber**



**Soft**

- Fashionable opaque hosiery for professional and everyday wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Cotton-lined gusset for pantyhose styles

### Fiber Content

- 79% Nylon , 21% Spandex  
(Content can vary slightly by style)
- Non-latex

### Best Suited for

- Management of early stages of chronic venous disease
- Pre- and post-vein treatment with minimal edema
- Patients who prefer fashionable options
- Pregnancy
- Long distance travel




**Silicone**  
band grip-top  
in calf



**Sensinno**  
band grip-top  
in thigh

20-30mmHg  
Closed Toe Thigh Medium Long in Mist

**842NMLW92**



30-40mmHg  
Open Toe Pantyhose Large Long in Nude

**843PML035**

**Sizing Chart**

	Small		Medium		Large	
Ankle Circumference	7"-8.5" (18-21.5cm)		8.5"-10" (21.5-25.5cm)		10"-11.5" (25.5-29cm)	
Calf Circumference	11"-15.5" (28-39cm)		12.5"-17.5" (31.5-45cm)		14"-20" (35.5-51cm)	
Calf Length	Short	Long	Short	Long	Short	Long
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
Thigh Circumference	18"-26" (45.5-66cm)		19"-28.5" (48-73cm)		20.5"-33" (52cm-84cm)	
Hip Circumference Guidelines	32.5"-44.5" (83-112.5cm)		34.5"-50" (87-127cm)		38.5"-59.5" (98-150.5cm)	
Leg Length	Short	Long	Short	Long	Short	Long
	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes Short leg 5.5-9.5, Long leg 8-12

Maternity hip circumference guidelines Small: <56" (142cm), Medium: <62" (157.5cm), Large: <70" (178cm)

**Styles and Colors**

	MIDNIGHT BLUE 09	NUDE 35	PECAN 40	TEAL 51*	MULBERRY 60*	BLUE IRIS 68*	CARAMEL 83*	ESPRESSO 89**
Calf (C)	●	●○	●	●	●	●	●	●
Calf w/Grip-Top (C/S)**								
Thigh w/Grip-Top (N) Sensinnov band grip-top	●	●○	●	●	●	●	●	●
Pantyhose (P)	●	●○	●	●	●	●	●	●
Maternity Pantyhose (M)*								

\* 15-20 & 20-30mmHg ONLY \*\* 20-30mmHg ONLY

	GRAPHITE 91	MIST 92*	BLACK 99
Calf (C)	●	●	●○
Calf w/Grip-Top (C/S)**			●
Thigh w/Grip-Top (N) Sensinnov band grip-top	●	●	●○
Pantyhose (P)	●	●	●○
Maternity Pantyhose (M)*			●

Product Code Example: **842CMSW68**

Series	mmHg Level	Style	Size	Gender/Open Toe	Color
84	2	C	MS	W	68



# sigvaris





## MOTION

### Push Yourself

Why choose **Motion**?

- Cushioned sole for added comfort
- Durable thermoregulating wool
- Moisture-wicking yarns

#### Motion Products

Previous	Series Number	New	Page
CUSHIONED COTTON	<b>360</b>	Cushioned Cotton	32
TRAVERSE SOCKS	<b>412</b>	High Tech	34
MERINO WOOL	<b>422</b>	Thermoregulating Wool	36
DAILY COMFORT	<b>850</b>	Comfort	38

360

## Motion CUSHIONED COTTON

for women & men

Motion Cushioned Cotton is soft and breathable with a padded sole for optimal comfort, and is ideal for sensitive skin.

### Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 58

### Features and Benefits



**Cushioned  
Sole**



**Cotton**



**Sensitive  
Skin**

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Flat, low-friction toe seam
- Diabetic friendly

### Fiber Content

- 64% Supima Cotton , 28% Nylon , 8% Spandex
- Non-latex

**SUPIMA**  
WORLD'S FINEST COTTONS

### Best Suited for

- Management of CVD symptoms
- People who stand for prolonged periods of time
- Patients with sensitive toes
- Work shoes and boots
- Casual shoes and sports shoes



20-30mmHg  
Closed Toe Calf Medium Long in Black

**362CMLM99**



20-30mmHg  
Closed Toe Calf Small Long in White

**362CSLW00**

**Sizing Chart**

	Small		Medium		Large		X Large	
Ankle Circumference	7"-8.5" (18-21.5cm)		8.5"-10" (21.5-25.5cm)		10"-11.5" (25.5-29cm)		11.5"-14" (29-36cm)	
Calf Circumference	11"-15.5" (28-39cm)		12.5"-17.5" (31.5-45cm)		14"-20" (35.5-51cm)		16"-24" (41-61cm)	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12    Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

**Styles and Colors**

Closed Toe = ●

WHITE 00	BLACK 99
-------------	-------------

Calf (C)

●	●
---	---

Product Code Example: **362CMLM00**

Series	mmHg Level	Style	Size	Gender	Color
36	2	C	ML	M	00

412

## Motion HIGH TECH

for women & men

Motion High Tech delivers optimal comfort where it counts with an anatomically-contoured shape and a cushioned sole.

### Compression Level

20-30mmHg

### Features and Benefits



**Moisture  
Wicking**



**Anatomically  
Contoured**



**Color**

- High-tech fibers for thermal control
- Cushioned footbed, toe, and Achilles tendon zones for optimal comfort
- Additional compression in calf improves circulation

### Fiber Content

- 64% Nylon, 17% Polyester, 19% Spandex
- Non-latex

### Best Suited for

- Management of symptoms related to chronic venous disease
- People who stand on their feet all day
- Walking and everyday activities
- Casual and occupational footwear



20-30mmHg  
Closed Toe Calf Small Medium in Limeade

**412CSM64**



20-30mmHg  
Closed Toe Calf Medium Large in Blood Orange

**412CML63**

**Sizing Chart**

	Small			Medium			Large			X Large		
Ankle Circumference	7"-8.5" (18-21.5cm)			8.5"-10" (21.5-25.5cm)			10"-11.5" (25.5-29cm)			11.5"-13" (29-33cm)		
Calf Circumference	11"-15.5" (28-39cm)			12.5"-17.5" (31.5-45cm)			14"-20" (35.5-51cm)			16"-23" (40.5-58cm)		
Women's Shoe Size	4.5-7.5	8-10.5	11-13.5	4.5-7.5	8-10.5	11-13.5	8-10.5	11-13.5	>14	8-10.5	11-13.5	>14
Men's Shoe Size	≤6.5	7-9.5	10-12.5	≤6.5	7-9.5	10-12.5	7-9.5	10-12.5	13-15.5	7-9.5	10-12.5	13-15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

**Styles and Colors**

	WHITE 00	BLUE 50	LIME 54	PINK 56	RED 59	STEEL BLUE 62	BLOOD ORANGE 63	LIMEADE 64	BLACK 99
Calf (C)	●	●	●	●	●	●	●	●	●

Product Code Example: **412CML64**

Series	mmHg Level	Style	Size	Color
41	2	C	ML	64

422

## Motion THERMOREGULATING WOOL

for women & men

Motion Thermoregulating Wool combines the benefits of natural Merino wool with cushioned soles, making it ideal for wherever the road takes you!

### Compression Level

20-30mmHg

Also available in 15-20mmHg, see page 74

### Features and Benefits



Wool



Cushioned Thermoregulating  
Sole



- Natural properties of wool make these socks odor-resistant
- Cushioned, padded soles help prevent blisters
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

### Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Non-latex

### Best Suited for

- Management of symptoms related to chronic venous disease
- People who stand on their feet all day
- Work boots, steel-toed shoes, and most occupational footwear



30-40mmHg  
Closed Toe Calf Large Long in Charcoal

**422CLL12**



30-40mmHg  
Closed Toe Calf Small Long in Olive

**422CSL31**



**Sizing Chart**

	Small			Medium			Large		
Ankle Circumference	7"-8.5" (18-21.5cm)			8.5"-10" (21.5-25.5cm)			10"-11.5" (25.5-29cm)		
Calf Circumference	11"-15.5" (28-39cm)			12.5"-17.5" (31.5-45cm)			14"-20" (35.5-51cm)		
Women's Shoe Size	4.5-7.5	8-10.5	11-13.5	4.5-7.5	8-10.5	11-13.5	8-10.5	11-13.5	>14
Men's Shoe Size	≤6.5	7-9.5	10-12.5	≤6.5	7-9.5	10-12.5	7-9.5	10-12.5	13-15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL

**Styles and Colors**

Closed Toe = ●



Calf (C)	●	●
----------	---	---

Product Code Example: **422CML12**

Series	mmHg Level	Style	Size	Color
42	2	C	ML	12

850

## Motion COMFORT

for women & men

Motion Comfort offers cushioned soles and moisture wicking to keep your legs feeling energized all day long.

### Compression Level

20-30mmHg

### Features and Benefits



**Cushioned  
Sole**



**Durable**



**Moisture  
Wicking**

- NILIT® Aquarius fibers provide moisture wicking to keep feet dry and comfortable
- Cushioned soles to absorb impact and help prevent blisters
- Smooth, flat toe seam

### Fiber Content

- 86% Nylon, 14% Spandex
- Non-latex

### Best Suited for

- Management of symptoms related to chronic venous disease
- Healthcare professionals
- Uniforms requiring black or solid color
- Individuals working in physically-demanding occupations



20-30mmHg  
Closed Toe Calf Large Long in Black

**852CLLM99**



20-30mmHg  
Closed Toe Calf Medium Long in Black

**852CMLW99**


**Sizing Chart**

	Small		Medium		Large	
Ankle Circumference	7"-8.5" (18-21.5cm)		8.5"-10" (21.5-25.5cm)		10"-11.5" (25.5-29cm)	
Calf Circumference	11"-15.5" (28-39cm)		12.5"-17.5" (31.5-45cm)		14"-20" (35.5-51cm)	
Calf Length	Short	Long	Short	Long	Short	Long
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12    Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

**Styles and Colors**

Closed Toe = ●



Calf (C) ●



Product Code Example: **852CMLM99**

Series	mmHg Level	Style	Size	Gender	Color
85	2	C	ML	M	99

# sigvaris



## SPECIALTY

### Take care of Yourself

Why choose **Specialty**?

Created for specific medical needs, including:

- Advanced venous conditions
- Wound care management
- Lymphatic support
- Diabetes needs

#### Specialty Products

Product	Series Number	Page
Natural Rubber	<b>500</b>	42
Secure	<b>550</b>	44
Secure Arm	<b>560</b>	46
Diabetic Compression Sock	<b>602</b>	48
Advance Armsleeves	<b>910</b>	50
Anti-Embolism	<b>930</b>	52

500

## Specialty

### NATURAL RUBBER

for women & men

Natural Rubber offers patients a durable, long-term containment solution for compression therapy with minimal fiber fatigue throughout the day.

#### Compression Level

30-40mmHg

40-50mmHg

50-60mmHg

#### Features and Benefits



**Natural Rubber**



**Durable**



**Wide Size Range**


- Ideal for patients with more serious indications
- Extended sizing for a custom-like fit

#### Fiber Content

- 70% Nylon, 30% Natural Latex Rubber (Calf & Thigh-high styles)
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibers (Pantyhose style)
- **Caution:** contains natural rubber latex that may cause allergic reaction

#### Best Suited for

- Severe varicosities
- Chronic edema post DVT
- Severe venous stasis
- Prevention of ulcer recurrence
- Early stages of lymphedema
- Active venous ulcers



**Silicone**  
band grip-top  
in thigh

30-40mmHg  
Open Toe Thigh Small Average Long in Beige

**503NS2077**



40-50mmHg  
Open Toe Calf Medium Full Long in Beige

**504CM4077**



**Sizing Chart**

	Small				Medium				Large				X Large			
<b>Ankle Circumference</b>	7.5"-9" (19-23cm)				9"-12" (23-30cm)				12"-15.5" (30-40cm)				15.5"-18" (40-46cm)			
<b>Calf Circumference</b>	Average		Full		Average		Full		Average		Full		Average		Full	
	11"-15" (28-38cm)		<17" (43cm)		12.5"-17" (32-43cm)		<19" (48cm)		14.5"-19" (36-48cm)		<22.5" (57cm)		18"-22.5" (45-57cm)		<27" (68cm)	
<b>Calf Circumference for 505C (Full only)</b>	505C: full calf only		12.5"-17" (32-43cm)		505C: full calf only		14.5"-19" (36-48cm)		505C: full calf only		18"-22.5" (45-57cm)		505C: not available in XL			
<b>Calf Length</b>	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)
<b>Thigh Circumference</b>	Average		Full		Average		Full		Average		Full		Custom sizes available. Please contact customer care for details.			
	16"-21" (41-54cm)		<24" (61cm)		19"-24" (48-61cm)		<28" (71cm)		23"-28" (58-71cm)		<32" (81cm)					
<b>Thigh Circumference for 505T (Full only)</b>	505T: full only		<24" (61cm)		505T: full only		<28" (71cm)		505T: full only		<32" (81cm)					
<b>Leg Length</b>	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long				
	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)				
<b>Hip Circumference Guidelines</b>	<48" (122cm)				<54" (137cm)				<61" (155cm)							
<b>Size</b>	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4

**Styles and Colors**

Open Toe = **O**



**Caution:** this product contains natural rubber latex which may cause allergic reactions. For all indications & contraindications, please see Contraindications & Cautions on page 95. Sigvaris accepts no liability for non-observance of contraindications or cautions.

<b>Calf (C)</b>	<b>O</b>	Calf only available in <b>S, M,</b> and <b>L</b> full circumference for <b>50-60mmHg</b>
<b>Thigh w/out Grip-Top (T)</b>	<b>O</b>	
<b>Thigh w/Grip-Top (N)</b> Silicone band grip-top	<b>O</b>	Not available in <b>50-60mmHg</b>
<b>Pantyhose (P)</b>	<b>O</b>	Pantyhose only available in <b>30-40mmHg</b> and only in full thigh circumference
<b>Thigh w/Waist Attachment (W)</b> Indicate leg <b>LEFT (/L)</b> or <b>RIGHT (/R)</b>	<b>O</b>	Thigh w/ waist attachment is available in average and full thigh circumferences in <b>30-40mmHg</b> and in full thigh circumferences only in <b>40-50mmHg</b>

Product Code Example: **504CS2077**

Series	mmHg Level	Style	Size	Gender/ Open Toe	Color
50	4	C	S2	O	77

Product Code Example: **503WL1077/L**

Series	mmHg Level	Style	Size	Gender/ Open Toe	Color	Leg
50	3	W	L1	O	77	/L

550

## Specialty

### SECURE

for women & men

Secure combines compression and containment with comfort and durability that's optimal for comfortably managing advanced venous edema, lymphedema, and wound care.

### Compression Level

20-30mmHg

30-40mmHg

40-50mmHg\*

### Features and Benefits



**Containment**



**Dotted Grip-top Band**



**Opaque**

- More graduated compression and containment than other traditional garments
- Durable double-covered in-lay yarns allow comfortable all-day wear

### Fiber Content

- 78% Nylon, 22% Spandex (Calf style)
- 76% Nylon, 24% Spandex (Thigh style)
- Non-latex

### Indications

- Ideal for patients with chronic venous insufficiency (CVI)
- Lymphedema / Lipedema
- Recurrent venous ulcers / Severe varicosities
- Deep vein thrombosis (DVT)
- Superficial phlebitis
- Post-thrombotic syndrome
- Orthostatic hypotension and post-phlebotic syndrome



Silicone band grip-top

30-40mmHg  
Closed Toe Calf Large Long in Black

**553CL2M99/S**



40-50mmHg  
Closed Toe Thigh Medium Full Long in Beige

**554NM4W77**

\* 40-50mmHg products made on demand.

**Sizing Chart**

	Medium				Large				X Large				XX Large			
<b>Ankle Circumference</b>	8.5"-10" (21.5-25.5cm)				10"-11.5" (25.5-29cm)				11.5"-13" (29-33cm)				13"-14.5" (33-37cm)			
<b>Calf Circumference</b>	Average		Full		Average		Full		Average		Full		Average		Full	
	12"-17.5" (31-45cm)		17.5"-23" (45-59cm)		13.5"-19" (35-49cm)		19"-25" (49-63cm)		16"-22" (41-56cm)		22"-28" (56-71cm)		17.5"-23" (45-59cm)		23"-29" (59-74cm)	
<b>Calf Length</b>	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)	≤16" (41cm)	>16" (41cm)
<b>Thigh Circumference</b>	Average		Full		Average		Full		Average		Full		Average		Full	
	20"-27" (51-69cm)		23.5"-31" (60-78cm)		24"-30" (61-77cm)		26.5"-34" (68-86cm)		27"-33.5" (69-85cm)		30"-37" (76-94cm)		29.5"-36.5" (75-93cm)		33"-40" (84-102cm)	
<b>Leg Length</b>	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)	≤28.5" (72cm)	>28.5" (72cm)
<b>Size</b>	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4	G1	G2	G3	G4

**Women's Shoe Sizes:** Short Leg 5.5-9.5, Long Leg 8-12    **Men's Shoe Sizes:** Short Leg 7-11, Long Leg 9.5-14

**Styles and Colors**

**Closed Toe = ●**  
**Open Toe = ○**



<b>Calf w/Grip-Top (C/S)</b> Silicone band grip-top	●○	●●
<b>Thigh w/Grip-Top (N)</b> Silicone band grip-top	●○	●●

**Product Code Example: 552CM2W77/S**

Series	mmHg Level	Style	Size	Gender/ Open Toe	Color	Grip-top
55	2	C	M2	W	77	/S

## Specialty SECURE ARMSLEEVE, GLOVE, AND GAUNTLET

for women & men

The Secure Armsleeve combines compression and containment with comfort and durability that's optimal for comfortably managing and helping to prevent lymphedema.

### Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

### Features and Benefits



**Containment**



**Dotted  
Grip-top Band**



**Colors**

- More graduated compression and containment than traditional garments
- Full-fingered Glove and Gauntlet available
- Provides UPF 50+ sun protection

### Fiber Content

- 78% Nylon, 22% Spandex (Armsleeve)
- 60% Nylon, 40% Spandex (Glove, Gauntlet)
- Non-latex

### Indications

- Ideal for patients with chronic edema
- Lymphedema
- Post-operative edema
- Post-mastectomy care



Silicone  
band grip-top

30-40mmHg  
Armsleeve Medium Short Average in Cocoa

**563AM84**



20-30mmHg  
Gauntlet Medium in Beige


**562GM78**

**Sizing Chart Arm**

	Small				Medium				Large				X Large				XX Large			
Wrist	5.5"-6.5" (14-16cm)				6.5"-7" (16-18cm)				7"-8" (18-20cm)				8"-9" (20-23cm)				9"-10" (23-26cm)			
Elbow	Average		Full		Average		Full		Average		Full		Average		Full		Average		Full	
	8"-10" 20-26cm		10"-13.5" 26-34cm		9"-11.5" 23-29cm		11.5"-14.5" 29-37cm		10"-12.5" 26-32cm		12.5"-15.5" 32-40cm		11.5"-13.5" 29-34cm		13.5"-16.5" 34-42cm		12.5"-14.5" 32-37cm		14.5"-17.5" 37-45cm	
Upper Arm	Average		Full		Average		Full		Average		Full		Average		Full		Average		Full	
	10"-13.5" 26-34cm		13.5"-15.5" 34-40cm		11.5"-14.5" 29-37cm		14.5"-17" 37-43cm		12.5"-15.5" 32-40cm		15.5"-18" 40-46cm		13.5"-16.5" 34-42cm		16.5"-19" 42-48cm		14.5"-17.5" 37-45cm		17.5"-20.5" 45-52cm	
Arm Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)
Size	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4	G1	G2	G3	G4

**Sizing Chart Hand**

	X Small		Small		Medium		Large		X Large		XX Large	
Palm Circumference	6"-6.5" (16-17cm)		7"-7.5" (18-19cm)		8"-8.5" (20-21cm)		8.5"-9" (22-23cm)		9.5"-10" (24-25cm)		10.5"-11" (26-28cm)	
Size	XS		S		M		L		XL		G	



**Armsleeve**

Product Code Example: **563AM278**

Series	mmHg Level	Style	Size	Color
56	3	A	M2	78

**Glove**

Product Code Example: **562FM78**

Series	mmHg Level	Style	Size	Color
56	2	F	M	78

**Gauntlet**

Product Code Example: **562GM78**

Series	mmHg Level	Style	Size	Color
56	2	G	M	78

**Styles and Colors**

Sleeve, Glove, and Gauntlet = ●



Armsleeve (A)	●	●	●	●
Glove (F)*	●	●	●	●
Gauntlet (G)*	●	●	●	●

\* Glove and Gauntlet available in 15-20mmHg & 20-30mmHg



602

## Specialty

### DIABETIC COMPRESSION SOCKS

for women & men

Diabetic Compression Socks are the only patented, clinically-proven graduated compression socks constructed for people with diabetes.



### Compression Level

18-25mmHg

Also available in 8-15mmHg, see page 78

### Features and Benefits



**Diabetes**



**Moisture Wicking**



**Cushioned Sole**

- Padded soles and invisible toe seam add comfort and protection
- Compression-free foot for optimum comfort
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- *drirelease*® yarns add moisture wicking properties

### Fiber Content

- 76% *drirelease* (65% Polyester, 11% Cotton), 19% Nylon, 5% Spandex
- Non-latex

### Best Suited for

- Individuals with diabetes and symptoms related to chronic venous disease and leg edema
- Individuals with peripheral neuropathy and small vessel disease



18-25mmHg  
Closed Toe Calf Medium Long in White

**602CMLM00**



18-25mmHg  
Closed Toe Calf Small Short in White

**602CSSW00**



**Sizing Chart**

	Small		Medium		Large		X Large	
<b>Ankle Circumference</b>	<b>7"-8.5"</b> (18-21.5cm)		<b>8.5"-10"</b> (21.5-25.5cm)		<b>10"-11.5"</b> (25.5-29cm)		<b>11.5"-14"</b> (29-36cm)	
<b>Calf Circumference</b>	<b>11"-15.5"</b> (28-39cm)		<b>12.5"-17.5"</b> (31.5-45cm)		<b>14"-20"</b> (35.5-51cm)		<b>16"-24"</b> (41-61cm)	
<b>Calf Length</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>	<b>Short</b>	<b>Long</b>
	<b>≤16"</b> (41cm)	<b>&gt;16"</b> (41cm)	<b>≤16"</b> (41cm)	<b>&gt;16"</b> (41cm)	<b>≤16"</b> (41cm)	<b>&gt;16"</b> (41cm)	<b>≤16"</b> (41cm)	<b>&gt;16"</b> (41cm)
<b>Size</b>	<b>SS</b>	<b>SL</b>	<b>MS</b>	<b>ML</b>	<b>LS</b>	<b>LL</b>	<b>XS*</b>	<b>XL*</b>

**Women's Shoe Sizes:** Short Leg 5.5-9.5, Long Leg 8-12    **Men's Shoe Sizes:** Short Leg 7-11, Long Leg 9.5-14

\*XS & XL Only Available In Men's Sizing

**Styles and Colors**

Closed Toe = ●



Calf (C)

**Product Code Example: 602CMLW00**

Series	mmHg Level	Style	Size	Gender/ Open Toe	Color
<b>60</b>	<b>2</b>	<b>C</b>	<b>ML</b>	<b>W</b>	<b>00</b>

## Specialty ADVANCE ARMSLEEVE

for women & men

The Advance Armsleeve combines the ultimate softness with silver ion properties and a Sensinnov grip-top band that's ideal for all-day comfort and effectiveness.

### Compression Level

15-20mmHg

20-30mmHg

30-40mmHg

### Features and Benefits



**Moisture  
Wicking**



**Sensinnov  
Grip-top  
Band**



**Breathable**


- Ideal for patients with arm lymphedema or post-mastectomy care
- Sensinnov top-band ensures armsleeve stays in place all day long
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable to provide comfortable, year-round wear

### Fiber Content

- 50% Nylon, 28% Rayon, 22% Spandex (Content can vary slightly by style)
- Non-latex

### Best Suited for

- Chronic edema post DVT
- Early stages of lymphedema
- Active venous ulcers



**Sensinnov  
band  
grip-top**

15-20mmHg  
Armsleeve Without Gauntlet  
Small Long in Beige

**911ASLO77/S**



**Sensinnov  
band  
grip-top**

20-30mmHg  
Armsleeve With Gauntlet  
Small Long in Beige

**912ASLG77/S**


**Sizing Chart**

	Small		Medium				Large				X Large			
Palm	6.5"-8" (17-21cm)		7"-8.5" (18-22cm)				7.5"-9" (19-23cm)				8"-9.5" (20-24cm)			
Wrist	5.5"-6.5" (14-16cm)		6.5"-7" (16-18cm)				7"-8" (18-20cm)				8"-9" (20-23cm)			
Upper Mid-Arm	Regular		Regular		Plus		Regular		Plus		Regular		Plus	
	10"-13.5" 26-34cm		11.5"-14.5" 29-37cm		14.5"-17" 37-43cm		12.5"-15.5" 32-40cm		15.5"-18" 40-46cm		13.5"-16.5" 34-42cm		16.5"-19" 42-48cm	
Arm Length	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long
	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)	≤17" (43cm)	>17" (43cm)
Size	SR	SL	MR	ML	MR/P	ML/P	LR	LL	LR/P	LL/P	XR	XL	XR/P	XL/P

**Styles and Colors**

Sleeve = ●

	
Armsleeve w/Gauntlet (G) Sensinnov band grip-top	●
Armsleeve w/out Gauntlet (O) Sensinnov band grip-top	●



**Note:** all armsleeves come with Grip-top (/S)  
 15-20mmHg not available with gauntlet  
 and not available in plus sizing

**Armsleeve without gauntlet**  
 Product Code Example: **912ASLO77/S**

Series	mmHg Level	Style	Size	W/out Gauntlet	Color	Grip-Top
91	2	A	SL	O	77	/S

**Armsleeve with gauntlet**  
 Product Code Example: **912ALRG77/PS**

Series	mmHg Level	Style	Size	W/ Gauntlet	Color	Plus Size	Grip-Top
91	2	A	LR	G	77	/P	S

930

## Specialty

### ANTI-EMBOLISM STOCKINGS

for women & men

Anti-Embolism Stockings are optimal for post-operative care and for combating Deep Vein Thrombosis (DVT) while on bed rest.

### Compression Level

18-23mmHg

### Features and Benefits



**Anatomical Foot Design**



**Breathable**



**Non-Ambulatory**

- Unique open-toe design allows for easy inspection of feet
- Anatomical foot design makes it easy to wear from bed to physical therapy
- Non-latex to prevent allergic reactions

### Fiber Content

- 71% Nylon, 29% Spandex  
(Content can vary slightly by style)
- Non-latex

### Best Suited for

- Individuals needing post-operative compression
- Bed rest
- Temporary use when transitioning from recovery to early rehabilitation



18-23mmHg  
Open Toe Thigh Small in White

**932NS000**



18-23mmHg  
Open Toe Calf Large in White

**932CLO00**

**Sizing Chart**

	<b>Small</b>	<b>Medium</b>	<b>Large</b>	<b>X Large</b>	<b>XX Large</b>
<b>Ankle Circumference</b>	<b>7"-8.5"</b> (18-22cm)	<b>8.5"-10"</b> (22-26cm)	<b>10"-11.5"</b> (26-30cm)	<b>11.5"-13"</b> (30-33cm)	<b>11.5"-13"</b> (30-33cm)
<b>Calf Circumference</b>	<b>11.5"-15"</b> (30-38cm)	<b>13.5"-16.5"</b> (34-42cm)	<b>15.5"-18.5"</b> (39-47cm)	<b>17.5"-21"</b> (45-53cm)	<b>19.5"-23"</b> (50-58cm)
<b>Thigh Circumference</b>	<b>18"-23.5"</b> (46-60cm)	<b>21"-26.5"</b> (53-68cm)	<b>23.5"-30"</b> (60-76cm)	<b>26.5"-33"</b> (68-84cm)	<b>30"-36"</b> (76-92cm)
<b>Size</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>X</b>	<b>2</b>

**Calf Length:** <18" (46cm) **Thigh Length:** <31.5" (80cm)

**Styles and Colors**

**Open Toe = O**



<b>Calf (C)</b>	<input type="radio"/>
<b>Thigh (N)</b>	<input type="radio"/>

**Product Code Example: 932CMO00**

Series	mmHg Level	Style	Size	Gender/ Open Toe	Color
<b>93</b>	<b>2</b>	<b>C</b>	<b>M</b>	<b>O</b>	<b>00</b>

# sigvaris





## WELL BEING

### Care for Yourself

Why choose **WELL BEING**?

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease
- Prevention of swelling and varicose veins
- Simplified sizing

#### WELL BEING & CONSUMER HEALTH Products

Product	Series Number	Page
CASUAL COTTON	<b>146 / 186</b>	56
CUSHIONED COTTON	<b>142 / 182</b>	58
SEA ISLAND COTTON	<b>151 / 191</b>	60
ALL-SEASON MERINO WOOL	<b>152 / 192</b>	62
MICROFIBER SHADES	<b>143 / 183</b>	64
SHEER FASHION	<b>120</b>	66
SOFT SILHOUETTE LEGGINGS	<b>170</b>	68
BUSINESS CASUAL	<b>189</b>	70
ATHLETIC RECOVERY SOCKS	<b>401</b>	72
PERFORMANCE SLEEVES	<b>412</b>	74
MERINO OUTDOOR SOCKS	<b>421</b>	76
EVERSOFT DIABETIC SOCK	<b>160</b>	78

146/186

# CASUAL COTTON

*for women & men*

## Compression Level

---

15-20mmHg

## Features and Benefits

---

- » Made with premium Supima Cotton for optimum softness and durability
- » Flat, low-friction toe seam to prevent irritation
- » Soft, stretchable knit-in band prevents pinching

## Fiber Content

---

- » 66% Supima Cotton, 26% Nylon, 8% Spandex
- » Non-latex

## Best Suited for

---

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease




**SUPIMA®**  
WORLD'S FINEST COTTONS

Sizing Chart

<b>146</b>	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	7"-9" (18-23cm)	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)
<b>186</b>	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)	8.5"-11.5" (21-29cm)
	Size	A	B	C

Style and Colors

Closed Toe = ●



Calf (C)	●	●	●	●	●
----------	---	---	---	---	---

Product Code Example: **146CB10**

Series	Style	Size	Color
<b>146</b>	<b>C</b>	<b>B</b>	<b>10</b>



**146**  
Women's Calf Style  
15-20mmHg  
Closed Toe Calf Size B in Khaki

**146CB30**



**186**  
Men's Calf Style  
15-20mmHg  
Closed Toe Calf Size C in Black

**186CC99**

# CUSHIONED COTTON

*for women & men*

## Compression Level

---

15–20mmHg

Also available in 20–30mmHg, see page 32

## Features and Benefits

---

- » Cushioned sole, heel, and toe area to protect feet
- » Flat, low-friction toe seam
- » Soft, stretchable knit-in band prevents pinching
- » Ideal for everyday wear and for sensitive skin
- » Suitable for patients with diabetes
- » No compression in heel and instep for ease of donning

## Fiber Content

---

- » 68% Supima Cotton, 25% Nylon, 7% Spandex
- » Non-latex

## Best Suited for

---

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease
- » Sensitive skin



**SUPIMA®**  
WORLD'S FINEST COTTONS

Sizing Chart

<b>142</b>	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	7"-9" (18-23cm)	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)
<b>182</b>	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)	8.5"-11.5" (21-29cm)
Size		A	B	C

Style and Colors

Closed Toe = ●



Product Code Example: **142CB00**

Series	Style	Size	Color
<b>142</b>	<b>C</b>	<b>B</b>	<b>00</b>



**142**  
Women's Calf Style  
15-20mmHg  
Closed Toe Calf Size B in White

**142CB00**



**182**  
Men's Calf Style  
15-20mmHg  
Closed Toe Calf Size C in Black

**182CC99**

# SEA ISLAND COTTON

*for women & men*

## Compression Level

---

15-20mmHg

Also available in 20-30mmHg, see page 14

## Features and Benefits

---

- » Knitted from the world's rarest cotton
- » Sea Island Cotton yarns are lightweight and have great strength, a silky luster, and an incredibly soft feel
- » Breathable and comfortable for sensitive skin
- » Ideal for travel and everyday wear

## Fiber Content

---

- » 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- » Non-latex

## Best Suited for

---

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease






Sizing Chart				
<b>151</b>	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	7"-9" (18-23cm)	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)
<b>191</b>	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)	8.5"-11.5" (21-29cm)
	Size	A	B	C

Style and Colors

Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---

Product Code Example: **151CA10**

Series	Style	Size	Color
<b>151</b>	<b>C</b>	<b>A</b>	<b>10</b>



**151**  
Women's Calf  
Style  
15-20mmHg  
Closed Toe Calf Size A in Navy

**151CA10**



**191**  
Men's Calf  
Style  
15-20mmHg  
Closed Toe Calf Size B in Black

**191CB99**

# ALL SEASON MERINO WOOL

*for women & men*

## Compression Level

---

15-20mmHg

Also available in 20-30mmHg, see page 16

## Features and Benefits

---

- » Soft, comfortable, and itch-free
- » Natural properties of wool make these socks odor resistant
- » Thermoregulating for year-round wear and keeps feet dry

## Fiber Content

---

- » 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- » Non-latex

## Best Suited for

---

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease




**ZURICH**  
COLLECTION

Sizing Chart				
<b>152</b>	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	7"-9" (18-23cm)	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)
<b>192</b>	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)	8.5"-11.5" (21-29cm)
	Size	A	B	C

Style and Colors

Closed Toe = ●



Calf (C)	●	●	●
----------	---	---	---

Product Code Example: **192CB11**

Series	Style	Size	Color
<b>192</b>	<b>C</b>	<b>B</b>	<b>11</b>



**152**  
Women's Calf  
Style

15-20mmHg  
Closed Toe Calf Size A in Black

**152CA99**



**192**  
Men's Calf  
Style

15-20mmHg  
Closed Toe Calf Size C in Navy

**192CC10**

# MICROFIBER SHADES

*for women & men*

## Compression Level

---

15–20mmHg

Also available in 20–30mmHg, see page 26

## Features and Benefits

---

- » Advanced microfiber technology for ultimate comfort
- » Glides on and off with ease
- » Excellent breathability and moisture management
- » Soft and durable sock for everyday wear

## Fiber Content

---

- » 87% Nylon, 13% Spandex
- » Non-latex

## Best Suited for

---

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease



**Sizing Chart**

<b>143</b>	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	7"-9" (18-23cm)	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)
<b>183</b>	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)	8.5"-11.5" (21-29cm)
Size		A	B	C

**Style and Colors**

Closed Toe = ●  
Actual color may vary



Calf (C)	●	●	●	●	●	●	●	●
----------	---	---	---	---	---	---	---	---

Closed Toe = ●  
Actual color may vary



Calf (C)	●	●	●	●	●	●
----------	---	---	---	---	---	---

**Product Code Example: 143CB49**

Series	Style	Size	Color
<b>143</b>	<b>C</b>	<b>B</b>	<b>49</b>



**143**  
Women's Calf  
Style

15-20mmHg  
Closed Toe Calf Size B in Purple Argyle

**143CB46**



**183**  
Men's Calf  
Style

15-20mmHg  
Closed Toe Calf Size C in Graphite Mini-Stripe

**183CC07**

# SHEER FASHION

*for women*

## Compression Level

15-20mmHg

## Features and Benefits

- » Fabric feels soft against the skin
- » Incredible sheer design in a variety of fashionable colors
- » Sheer toe and heel area
- » Comfortable stretch is excellent for everyday wear
- » Cotton-lined gusset on pantyhose styles

## Fiber Content

- » 79% Nylon, 21% Spandex
- » Content can vary slightly by style
- » Non-latex

## Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic vein disease



15-20mmHg  
Closed Toe Calf Size B in Mulberry

**120CB60**



15-20mmHg  
Closed Toe Pantyhose Size C in Black

**120PC99**



**Pantyhose and Maternity Sizing Chart**

Ft/Cm	Lbs/Kg	100/45	105/48	110/50	115/52	120/54	125/57	130/59	135/61	140/64	145/66	150/68	155/70	160/73	165/75	170/77	175/80	180/82	185/84	190/86	195/89	200/91	205/93	210/95	215/98	220/100	225/102	230/105
4'9" / 145cm																												
4'10" / 147cm			<b>SIZE A</b>																									
4'11" / 150cm			ANKLE CIRCUMFERENCE																									
5'0" / 152cm			7.5"-9" (19-23cm)																									
5'1" / 155cm																												
5'2" / 157cm																												
5'3" / 160cm																												
5'4" / 162cm			<b>SIZE B</b>																									
5'5" / 165cm			ANKLE CIRCUMFERENCE																									
5'6" / 167cm			7.5"-9" (19-23cm)																									
5'7" / 170cm																												
5'8" / 173cm																												
5'9" / 175cm																												
5'10" / 178cm																												
5'11" / 180cm																												
6' / 183cm																												

**Note:** please use pre-pregnancy weight for maternity pantyhose. \* Maternity pantyhose not available in size F.

**Calf / Thigh Sizing Chart**

<b>120</b>	Women's Shoe	5-7	7.5-9.5	10-12
	Ankle Circumference	7.5"-9" (19-23cm)	8"-10" (20.5-25cm)	8.5"-10.5" (22-27cm)
	Size	A	B	C

**Weight Guidelines For Thigh Style:** Size A: <135Lbs., Size B: <180Lbs., Size C: <230Lbs.

**Styles and Colors**

Closed Toe = ●  
Open Toe = ○

<b>CHARCOAL</b> 12	<b>TAUPE</b> 29	<b>NATURAL</b> 33	<b>SUNTAN</b> 36	<b>MULBERRY</b> 60	<b>CAFÉ</b> 73	<b>BLACK</b> 99
-----------------------	--------------------	----------------------	---------------------	-----------------------	-------------------	--------------------

Calf (C)	●	●●	●●	●●	●	●	●
Thigh (N)	●	●	●	●		●	●
Pantyhose (P)	●	●	●	●		●	●
Maternity Pantyhose (M)			●				●

**Product Code Example: 120PE73**

Series	Style	Size	Color
120	P	E	73

**Product Code Example: 120CCO36**

Series	Style	Size	Open Toe	Color
120	C	C	O	36

# SOFT SILHOUETTE LEGGINGS

*for women*

## Compression Level

---

15-20mmHg

## Features and Benefits

---

- » Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts
- » Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort
- » Advanced soft stretch microfiber fabric shapes legs and smooths skin

## Fiber Content

---

- » 83% Nylon, 17% Spandex
- » Non-latex

## Best Suited for

---

- » Tired, heavy, restless, aching legs,
- » Occupational leg health risks, family history of chronic vein disease



### Sizing Chart

Ft	Lbs																											
	100	105	110	115	120	125	130	135	140	145	150	155	160	165	170	175	180	185	190	195	200	205	210	215	220	225	230	
4'10"																												
4'11"																												
5'0"																												
5'1"																												
5'2"																												
5'3"																												
5'4"																												
5'5"																												
5'6"																												
5'7"																												
5'8"																												
5'9"																												
5'10"																												
5'11"																												
6'																												

**Guidelines for Ankle Circumference**  
 Size A: 7"-9", Size B & C: 8"-10", Size D: 8.5"-11.5"

### Style and Colors

Open Foot = <input type="radio"/>	<b>BLACK</b> 99
Leggings (L) <input checked="" type="radio"/>	<input type="radio"/>

**Product Code Example: 170LB99**

Series	Style	Size	Color
<b>170</b>	<b>L</b>	<b>B</b>	<b>99</b>



15-20mmHg Legging Size B in Black

**170LB99**



**Non-binding ankle cuff** stays comfortably in place for all day wear

# BUSINESS CASUAL

*for men*

## Compression Level

---

15-20mmHg

## Features and Benefits

---

- » Stretchable toe fabric and heel design increase comfort and durability
- » Classic ribbed look for everyday wear
- » Wide, non-binding top band for comfortable fit
- » Soft, microfiber inlay yarns

## Fiber Content

---

- » 88% Nylon, 12% Spandex
- » Non-latex

## Best Suited for

---

- » Tired, heavy, restless, aching legs,
- » Travelers, occupational leg health risks, family history of chronic vein disease



Sizing Chart				
189	Men's Shoe	6.5-8.5	9-11	11.5-14
	Ankle Circumference	7.5"-10" (19-25cm)	8"-10.5" (20-27cm)	8.5"-11.5" (21-29cm)
	Size	A	B	C

Style and Colors

Closed Toe = ●

	NAVY 10	BROWN 11	CHARCOAL 12	BLACK 99
Calf (C)	●	●	●	●

Product Code Example: 189CB12

Series	Style	Size	Color
189	C	B	12



15-20mmHg  
Closed Toe Calf Size B in Charcoal  
**189CB12**



15-20mmHg  
Closed Toe Calf Size C in Brown  
**189CC11**

# ATHLETIC RECOVERY SOCKS

for women & men

## Compression Level

---

15-20mmHg

## Features and Benefits

---

- » Helps reduce lactic acid buildup
- » Comfortable and easy to put on and take off
- » *drirelease* yarns provides moisture wicking and odor control
- » No compression in heel and instep for ease of donning

## Fiber Content

---

- » 67% *drirelease* (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- » Non-latex

## Best Suited for

---

- » Tired, heavy, restless, aching legs
- » Post-workout recovery
- » Mild athletic activity, walking





Sizing Chart

401	Women's Shoe	<7.5	7.5-9.5	10-12	12.5-15	≥15.5
	Men's Shoe	<6.5	6.5-8.5	9-11	11.5-14	14.5-17
	Ankle Circumference	7"-9.5" (18-24cm)	7.5"-10" (19-26cm)	8"-10.5" (20-27cm)	8.5"-11.5" (21-29cm)	8.5"-11.5" (21-29cm)
	Size	S	M	L	X	2

Style and Colors

Closed Toe = ●



Calf (C)



Product Code Example: 401C200

Series	Style	Size	Color
401	C	2	00



15-20mmHg  
Closed Toe Calf Size XX-Large in Black

**401C299**



15-20mmHg  
Closed Toe Calf Size Medium in White

**401CM00**

# PERFORMANCE SLEEVES

*for women & men*

## Compression Level

---

20–30mmHg

## Features and Benefits

---

- » Fabric resists deterioration in salt, chlorine, or fresh water
- » Stabilizes and supports active muscles
- » Quick drying
- » High-visibility reflective transfer

## Fiber Content

---

- » 75% Nylon, 25% Spandex
- » Non-latex

## Best Suited for

---

- » Running, triathlons, endurance sports, and gym workouts
- » Team sports
- » Occupational wear



Sizing Chart

<b>412</b>	Ankle Circumference	7"–8.5" (18–21.5cm)	8.5"–10" (21.5–25.5cm)	10"–11.5" (25.5–29cm)	11.5"–13" (29–33cm)
	Calf Circumference	11"–15.5" (28–39cm)	12.5"–17.5" (31.5–45cm)	14"–20" (35.5–51cm)	16"–23" (40.5–58cm)
	Size	S	M	L	X

Style and Colors

Calf Sleeve = ●	WHITE 00	BLUE 50	LIME 54	PINK 56*	PURPLE 58*	RED 59	STEEL BLUE 62	BLOOD ORANGE 63	LIMEADE 64	BLACK 99
Calf Sleeve (V)	●	●	●	●	●	●	●	●	●	●

\* Pink and purple not available in X-Large

Product Code Example: **412VL63**

Series	Style	Size	Color
<b>412</b>	<b>V</b>	<b>L</b>	<b>63</b>



20–30mmHg  
Calf Sleeve Size Large in Blood Orange

**412VL63**



20–30mmHg  
Calf Sleeve Size Medium in Steel Blue

**412VM62**

# MERINO OUTDOOR SOCKS

*for women & men*

## Compression Level

---

15–20mmHg

Also available in 20–30mmHg, see page 36

## Features and Benefits

---

- » Soft, comfortable, and itch-free
- » Natural properties of wool make these socks odor-resistant
- » Mesh flex zones for comfort and ventilation
- » Cushioned, padded soles protect against blisters
- » No compression in heel and instep for ease of donning
- » Thermoregulating for year-round comfort

## Fiber Content

---

- » 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- » Non-latex

## Best Suited for

---

- » Tired, heavy, restless, aching legs
- » Hiking, fishing, and hunting
- » Post-workout recovery
- » Outdoor sports and walking



Sizing Chart

421	Women's Shoe	<7.5	7.5-9.5	10-12	12.5-15
	Men's Shoe	<6.5	6.5-8.5	9-11	11.5-14
	Ankle Circumference	7"-9.5" (18-24cm)	7.5"-10" (19-26cm)	8"-10.5" (20-27cm)	8.5"-11.5" (21-29cm)
	Size	S	M	L	X

Style and Colors

Closed Toe = ●



Calf (C)	●	●
----------	---	---

Product Code Example: 421CS31

Series	Style	Size	Color
421	C	S	31



15-20mmHg  
Closed Toe Calf Size Large in Olive

**421CL31**



15-20mmHg  
Closed Toe Calf Size Medium in Charcoal

**421CM12**

# EVERSOFT DIABETIC SOCKS

*for women & men*

## Compression Level

---

8-15mmHg

Also available in 18-25mmHg, see page 48

## Features and Benefits

---

- » Thick, padded soles, heels, and toes for comfort and protection
- » Flat, low friction toe seam
- » Soft, stretchable knit-in band prevents pinching
- » Compression from ankle to knee, compression-free foot

## Fiber Content

---

- » 67% *drirelease* (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- » Non-latex

## Best Suited for

---

- » Individuals with diabetes
- » Individuals with limited dexterity or lower tolerance for compression
- » Individuals with sensitive skin





Sizing Chart

160	Women's Shoe	<7.5	7.5-9.5	10-12	12.5-15
	Men's Shoe	<6.5	6.5-8.5	9-11	11.5-14
	Ankle Circumference	7"-9.5" (18-24cm)	7.5"-10" (19-26cm)	8"-10.5" (20-27cm)	8.5"-11.5" (21-29cm)
	Size	S	M	L	X

Style and Colors

Closed Toe = ●



Calf (C)



Product Code Example: 160CS00

Series	Style	Size	Color
160	C	S	00



8-15mmHg  
Closed Toe Calf Size Small in White

**160CS00**



8-15mmHg  
Closed Toe Calf Size Large in White

**160CL00**

# sigvaris



# Medically Complex Edema

## MEDICALLY COMPLEX EDEMA

The SIGVARIS GROUP line of products for medically complex edema (MCE) are for the treatment and management of conditions impacting the lymphatic and venous circulatory systems.

The product portfolio contains ready-to-wear garments for the active treatment and maintenance of lymphedema and for patients managing chronic venous diseases in the upper and lower extremities.

A sampling of these products is included in the following pages.

For more details and to request a complete MCE catalog, please ask your local SIGVARIS GROUP Territory Manager or email [sigvariscares@sigvaris.com](mailto:sigvariscares@sigvaris.com).

### Medically Complex Edema Products

Product	Page
Compreflex Standard Calf & Foot	82
Compreflex Standard Calf	82
Compreflex Transition Calf	84
Compreflex Complete Calf	84
Coolflex Standard Calf	86
Compreflex Strap Extender	86
Assorted Liners	87
Assorted Foot Garments	88

## Compreflex® Standard CALF & FOOT

20-50mmHg



### When to Use

- Adjustable ranges of resting compression is desired (20-50mmHg)
- Non-attached foot piece is desired
- No compression in the liner is desired
- Patient requires a calf and foot garment

### Conditions

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lipedema



**Compreboot Standard**

30-40mmHg  
(see page 88 for details)



**Basic Liner**

Non-Compressive  
(see page 87 for details)



### Accutab® Devices

Accutab Devices  
(visit website for details)

## Compreflex Standard CALF

20-50mmHg



### When to Use

- Varying ranges of adjustable resting compression is desired (20-50mmHg)
- No compression in the liner is desired
- No Compreboot needed
- Patient requires compression in the calf only or a different foot option (size or boot)
- Ability to wear standard footwear over liner is desired

### Conditions

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lipedema



**Basic Liner**

Non-Compressive  
(see page 87 for details)



### Accutab® Devices

Accutab Devices  
(visit website for details)




30–40mmHg Compreboot for added foot and ankle compression (see page 88 for details)

## Compreflex Standard Calf & Foot Sizing Chart

	Small	Medium	Large	X Large	XX Large	
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm	
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm	
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm	
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm	
Black	Regular	1401-BKR	1402-BKR	1403-BKR	1404-BKR	1405-BKR
	Tall	1401-BKT	1402-BKT	1403-BKT	1404-BKT	1405-BKT
Beige	Regular	1411-BKR	1412-BKR	1413-BKR	1414-BKR	1415-BKR
	Tall	1411-BKT	1412-BKT	1413-BKT	1414-BKT	1415-BKT

Tall garments ship with Regular length Compreboot. If Long Compreboot is required, purchase Compreflex Standard Calf and order desired Compreboot separately.

## REQUIRED MEASUREMENTS



**K:** Measure length on medial side from heel to 1st metatarsal head

**G Length**  
Regular: 30–33.9cm  
Tall: 34–42cm

## Compreflex Standard – Calf Sizing Chart (A–C1)

	Small	Medium	Large	X Large	XX Large	
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm	
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm	
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm	
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm	
Black	Regular	1401-NF-BKR	1402-NF-BKR	1403-NF-BKR	1404-NF-BKR	1405-NF-BKR
	Tall	1401-NF-BKT	1402-NF-BKT	1403-NF-BKT	1404-NF-BKT	1405-NF-BKT
Beige	Regular	1411-NF-BKR	1412-NF-BKR	1413-NF-BKR	1414-NF-BKR	1415-NF-BKR
	Tall	1411-NF-BKT	1412-NF-BKT	1413-NF-BKT	1414-NF-BKT	1415-NF-BKT



For a PDF of the product ordering form, scan the QR code or visit [bit.ly/SIGmccOrderforms](https://bit.ly/SIGmccOrderforms)



## Compreflex Transition

### CALF

20-50mmHg



#### When to Use

- Medicare approved for ulcers under code **A6545-AW**
- Patient requires a compressive liner; the Transition Liner provides 10-15mmHg in the foot and ankle
- Mild foot compression in liner is desired
- Ability to wear standard footwear over liner is desired

#### Conditions

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lipedema



**Transition Liner**

10-15 mmHg of mild compression in the foot and ankle.  
(see page 87 for details)



#### Accutab Devices

Accutab Devices  
(visit website for details)

## Compreflex Complete

### CALF

20-50mmHg



#### When to Use

- Patient requires a compressive liner; the Complete Liner provides 10-15mmHg in the foot, ankle, and calf
- Mild foot and calf compression in liner is desired
- Ability to wear standard footwear over liner is desired

#### Conditions

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lipedema



**Complete Liner**

Mild compression in foot, ankle, and calf  
(see page 87 for details)




#### Accutab Devices

Accutab Devices  
(visit website for details)



## REQUIRED MEASUREMENTS



### Compreflex Transition – Calf Sizing Chart (A-C1)

	Small	Medium	Large	X Large	XX Large	
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm	
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm	
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm	
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm	
Black	Regular	1401-UC-BKR	1402-UC-BKR	1403-UC-BKR	1404-UC-BKR	1405-UC-BKR
	Tall	1401-UC-BKT	1402-UC-BKT	1403-UC-BKT	1404-UC-BKT	1405-UC-BKT
Beige	Regular	1411-UC-BKR	1412-UC-BKR	1413-UC-BKR	1414-UC-BKR	1415-UC-BKR
	Tall	1411-UC-BKT	1412-UC-BKT	1413-UC-BKT	1414-UC-BKT	1415-UC-BKT



**Compreflex Complete Calf**  
20-50mmHg  
Medium Regular in Beige  
**1412-FC-BKR**

### Compreflex Complete – Calf Sizing Chart (A-C1)

	Small	Medium	Large	X Large	XX Large	
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm	
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm	
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm	
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm	
Black	Regular	1401-FC-BKR	1402-FC-BKR	1403-FC-BKR	1404-FC-BKR	1405-FC-BKR
	Tall	1401-FC-BKT	1402-FC-BKT	1403-FC-BKT	1404-FC-BKT	1405-FC-BKT
Beige	Regular	1411-FC-BKR	1412-FC-BKR	1413-FC-BKR	1414-FC-BKR	1415-FC-BKR
	Tall	1411-FC-BKT	1412-FC-BKT	1413-FC-BKT	1414-FC-BKT	1415-FC-BKT



For a PDF of the product ordering form, scan the QR code or visit [bit.ly/SIGmccOrderforms](https://bit.ly/SIGmccOrderforms)

## Coolflex™ Standard CALF

20–40mmHg




### When to Use

- Varying ranges of adjustable resting compression is desired (20-40mmHg)
- No compression in the liner is desired
- Simplified one-handed donning is desired
- More breathable fabric is desired
- Desire to select separate foot piece

### Conditions

- Lymphedema
- Edema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lipedema



#### One-Handed Donning

Spacer Fabric TR is designed to be highly breathable for cool-wearing comfort

Durable, metal hook closures for maximum strength and durability



#### Basic Liner

Non-Compressive  
(see page 87 for details)



### Coolflex Standard Calf Sizing Chart (A–C1)

SIZE	Small	Small Max	Medium/ Large	Medium / Large Max	X Large / XX Large
C1	28–39cm	35–50cm	34–50cm	45–64cm	45–64cm
C	28–39cm	35–50cm	34–50cm	45–64cm	45–64cm
B	23–33cm	27–38cm	29–43cm	35–49cm	37–54cm
A	16–26cm	16–26cm	23–33cm	23–33cm	30–42cm

#### Left

Regular	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
Tall	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L

#### Right

Regular	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
Tall	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R

## Compreflex Strap Extender

- Provides up to 10cm of additional circumferential size range to Compreflex calf wraps
- Detaches easily
- Potential inexpensive solution for patients outside of standard size ranges of Compreflex calf wraps

	1 Unit	10 Unit Pack
Black	1401-SE1	1401-SEB
Beige	1411-SE1	1411-SEB



## Basic Liner

**Non-compressive**

### Features

- Non-compressive liner for the protection of the skin
- Easy to apply knee-high liners
- *drive/lease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Heathered Silver color only

Size: One size fits all	
Pair	10 Pair Pack
1700 - NK1	1700 - NKC
Included With Compreflex Standard Calf & Foot Compreflex Standard Calf	



## Transition Liner

**10-15mmHg**

### Features

- Mild compression in foot and ankle only
- *drive/lease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige

**Transition Liner Sizing Chart** (Included with Compreflex Transition)

SIZE	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe	Pair (Black)	Pair (Beige)
Small	19-26cm	29-39cm	<7½	<6½	1601-CK	1611-CK
Medium	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK	1612-CK
Large	26-36cm	40-50cm	10-12	9-11	1603-CK	1613-CK
X Large	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK	1614-CK
XX Large	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK	1615-CK



## Complete Liner

**10-15mmHg**

### Features

- Mild compression in foot and ankle with graduated compression in calf
- Comfortable knee-high liners
- *drive/lease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige

**Transition Liner Sizing Chart** (Included with Compreflex Complete)

SIZE	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe	Pair (Black)	Pair (Beige)
Small	19-26cm	29-39cm	<7½	<6½	1601-CK-F1	1611-CK-F1
Medium	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK-F1	1612-CK-F1
Large	26-36cm	40-50cm	10-12	9-11	1603-CK-F1	1613-CK-F1
X Large	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK-F1	1614-CK-F1
XX Large	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK-F1	1615-CK-F1



## Compreboot Lite

Variable compression



### Features

- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design (may be trimmed to fit)
- Ideal for mild edema/lymphedema
- Breathable stretch fabric conforms to ankle
- Provides mild, therapeutic compression
- Available in Black only



## Compreboot Standard

Variable compression



### Features

- 4-way stretch for firm compression around ankle
- Included with Compreflex Standard Calf & Foot
- Ideal for mild-moderate edema/lymphedema
- Features neoprene material with hook and loop tabs
- Open heel and toe
- Provides moderate, therapeutic compression
- Available in Beige and Black



## Compreboot Plus

Variable compression




### Features

- Chipped foam lining
- Ideal for chronic/acute edema/lymphedema
- Removable gripper patches provide traction
- Straps are removable for custom positioning or replacement
- Provides moderate, therapeutic compression
- Available in Black only



## REQUIRED MEASUREMENTS



## Measuring Instructions

**I:** Measure circumference of ankle bend and heel.

**J:** Measure circumference across metatarsal heads.

**K:** Measure length on the medial side of the foot from the heel to the 1st metatarsal head.

### Compreboot Lite and Compreboot Standard Chart

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–36cm	28–36cm	39cm max	39cm max	44cm max	44cm max
J	22–26cm	22–26cm	30cm max	30cm max	33cm max	33cm max
K	14–18cm	19–23cm	18–20cm	23–25cm	20–22cm	25–27cm

### Compreboot Lite Product Codes

		Small/ Medium/ Large	X Large XX Large
Black	Regular	1803-BLR	1805-BLR
	Long	1803-BLL	1805-BLL

### Compreboot Standard Product Codes

		Small	Medium Large	X Large XX Large
Black	Regular	1801-BTR	1803-BTR	1805-BTR
	Long	1801-BTL	1803-BTL	1805-BTL
Beige	Regular	1811-BTR	1813-BTR	1815-BTR
	Long	1811-BTL	1813-BTL	1815-BTL

### Compreboot Plus Chart

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–36cm	28–34cm	39cm max	39cm max	44cm max	44cm max
J	22–26cm	22–25cm	30cm max	28cm max	33cm max	33cm max
K	14–18cm	19–23cm	18–20cm	23–25cm	20–22cm	25–27cm

### Compreboot Plus Product Codes

		Small	Medium/ Large	X Large XX Large
Black	Regular	1801-BPR	1803-BPR	1805-BPR
	Long	1801-BPL	1803-BPL	1805-BPL



# sigvaris





## ACCESSORIES

### Help Yourself

Why choose **Accessories**?

- Aid in putting on and removing stockings
- Stocking cleaning to maximize usage
- Adhesives and extra liners for added comfort

#### Accessories

Product	Page
Gloves	92
Sim-Slide	92
Magnide On/Off	92
SIMON	93
Washing Solution	93
HOLD UP	93
Cotton Underliner	93
Foot Slip	94
Sigvaris Doff N' Donner	94
Sigvaris Doff N' Donner Cone	94

## GLOVES

**Rubber Gloves** – Ridged  
The practical and protective way to apply and remove your compression stockings.



Package: 12 per Pack

Size	Product Code
Small	<b>591R400S</b>
Medium	<b>591R400M</b>
Large	<b>591R400L</b>
X-Large	<b>591R400X</b>

Package: Single Pair

Size	Product Code
Small	<b>591RPRS</b>
Medium	<b>591RPRM</b>
Large	<b>591RPRL</b>
X-Large	<b>591RPRX</b>

**Latex-Free Gloves**  
Gripping donning gloves recommended for use in a non-latex environment.



Package: 12 per Pack

Size	Product Code
Small	<b>592R400S</b>
Medium	<b>592R400M</b>
Large	<b>592R400L</b>
X-Large	<b>592R400X</b>

Package: Single Pair

Size	Product Code
Small	<b>592RPRS</b>
Medium	<b>592RPRM</b>
Large	<b>592RPRL</b>
X-Large	<b>592RPRX</b>

Glove Sizing Chart			
Measure Circumference of the Hand at the Palm			
6.5"–7" (16.5–18cm)	7.5"–8" (19–20.5cm)	8.5"–9" (21.5–23cm)	9.5"–10" (24–25.5cm)
Small	Medium	Large	X-Large



For more **SIGVARIS** videos, including donning tips and product demonstrations, visit our YouTube channel: [youtube.com/c/sigvarisgroupusa](https://youtube.com/c/sigvarisgroupusa)

**Caution:** Rubber Gloves contains natural rubber latex that may cause allergic reaction

## DONNING AID



**Sim-Slide®**



**Video: How to Use the Sim-Slide**  
To watch the video go to: [bit.ly/SIGssyt](https://bit.ly/SIGssyt)

**Sim-Slide**

Size	Product Code
Small	<b>587SLIDESM</b>
Medium	<b>587SLIDEMD</b>
Large	<b>587SLIDELG</b>
X-Large	<b>587SLIDEXL</b>

**Sizing Chart**  
For use with Open toe styles

**NEW SIZING**

Women's Shoe Size	≤4–6	6.5–8	8.5–9.5	≥10
Men's Shoe Size	3–4.5	5–6.5	7–8	≥8.5
Size	S	M	L	XL



**Magnide On/Off®**



**Video: How to Use the Magnide On/Off**  
To watch the video go to: [bit.ly/SIGmofyt](https://bit.ly/SIGmofyt)

**Magnide On/Off**

Size	Product Code
Medium	<b>587MAGMD</b>
Large	<b>587MAGLG</b>
X-Large	<b>587MAGXL</b>

**Sizing Chart**  
For use with Closed toe styles

Women's Shoe Size	≤8.5	9–12	≥12.5
Men's Shoe Size	≤7.5	8–11	≥11.5
Size	M	L	XL*

\*XL suitable for extra-wide feet

# ACCESSORIES



## SIMON

A donning and doffing aid so versatile that it can easily be stowed to accompany you wherever you go.



*Video: How to Use the SIMON*

To watch the video go to: [bit.ly/SIGSIMONC](http://bit.ly/SIGSIMONC)

Size	Calf Circumference	Product code
X-Small	25–29cm	<b>589SXS</b>
Small	30–38cm	<b>589SSM</b>
Medium	38–46cm	<b>589SMD</b>
Large	46–54cm	<b>589SLG</b>



## Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for fine washables!

### Package Options

Contains	Product Code
Six 16 Oz. Bottles (473ml)	<b>586W900</b>
One 16 oz. Bottle (473ml)	<b>586W950</b>
Twelve 2 oz. Bottles (60ml)	<b>586W975</b>
One 2 oz. Bottle (60ml)	<b>586W985</b>



## Hold Up

Roll-on water-soluble liquid adhesive to keep stockings in place.

### Package Options

Contains	Product code
Twelve 2 oz. Bottles (60ml)	<b>5811300</b>
One 2 oz. Bottle (60ml)	<b>5811350</b>



## Cotton Underliner

10–15mmHg stocking. Can be used as an under layer to protect fragile skin, secure wound dressings, add 10–15mmHg of compression, or create an easy glide surface to aide in putting on and taking off compression stockings. For sizing, please refer to the Natural Rubber Calf sizing on page 42.

### Package: 2 per pack

Size	Product Code
Small/Short	<b>588USS</b>
Small/Long	<b>588USL</b>
Medium / Short	<b>588UMS</b>
Medium / Long	<b>588UML</b>
Large-XL / Short	<b>588ULS</b>
Large-XL / Long	<b>588ULL</b>



**Foot Slip**

Comfortable applicator for use with open toe stockings.

Package: 6 per pack

Size	Product Code
One Size	<b>580F200</b>



**SIGVARIS Doff N' Donner**

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use.

Size	Product Code
One Size	<b>583D100</b>



**SIGVARIS Doff N' Donner Cone**

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner. The bottom of the base firmly attaches to most clean, flat surfaces via suction.



*Video: How to Use the SIGVARIS Doff N' Donner*  
To watch the video go to:  
**bit.ly/SIGdndyt**

Size	Product Code
One Size	<b>583D200</b>



## SYMPTOMS & INDICATIONS

Symptoms / Indications	Compression Level at the Ankle
<ul style="list-style-type: none"> <li>Tired, achy legs</li> <li>Occupational or evening edema</li> <li>Leg discomfort from long hours of standing or sitting</li> <li>Prevent varicosities &amp; swelling during pregnancy</li> <li>Preventative for legs predisposed to risk</li> <li>Reduce swelling during travel</li> <li>Prevent DVT during travel</li> </ul>	<b>15–20mmHg</b>
<ul style="list-style-type: none"> <li>Diabetics who experience edema (see Cautions below)</li> </ul>	<b>18–25mmHg</b>
<ul style="list-style-type: none"> <li>Heavy, fatigued, aching legs</li> <li>Mild edema in lower extremities</li> <li>Mild varicosities with minimal edema</li> <li>Mild varicosities during pregnancy</li> <li>Lymphedema upper extremities</li> <li>Post-procedure of small veins</li> <li>Prevention or management of DVT</li> </ul>	<b>20–30mmHg</b>
<ul style="list-style-type: none"> <li>Moderate varicosities with mild to moderate edema (also during pregnancy)</li> <li>Post-procedure of larger veins to maintain treatment success</li> <li>CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations)</li> <li>Management and treatment of DVT or Post-Thrombotic Syndrome</li> <li>Superficial phlebitis</li> <li>Lymphedema after decongestant therapy to maintain reduction</li> <li>Orthostatic / postural hypotension</li> <li>Post-fracture, post-traumatic edema</li> </ul> <p><i>Note: A prescription is recommended for stockings with the indications of 30mmHg or higher</i></p>	<b>30–40mmHg</b>
<ul style="list-style-type: none"> <li>Lymphedema after decongestant therapy to maintain reduction</li> <li>Severe Post-Thrombotic Syndrome</li> <li>Severe varicosities and/or edema</li> </ul>	<b>40mmHg or higher</b>
<ul style="list-style-type: none"> <li>Severe arterial insufficiency, intermittent claudication, ischemia</li> <li>Uncontrolled congestive heart failure</li> </ul>	<b>Contraindications</b>
<ul style="list-style-type: none"> <li>Skin sensitivities or allergies</li> <li>Advanced neuropathy with sensory loss in leg or foot</li> <li>Diabetes with arterial compromise</li> <li>Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician</li> <li>Eczema, hypodermatitis, dermatosclerosis, dermatitis</li> </ul>	<b>Cautions</b>
<p>The therapeutic effect of graduated compression stockings is based on the compression level. With the proper CEAP classification, the correct compression level at the ankle can be chosen. For more severe symptoms, consult with a vascular specialist.</p>	
<p><i>No liability accepted for non-observance of contraindications and cautions.</i>            Fronek HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.</p>	

Term	Definition	Sigvaris Usage
CHRONIC VENOUS DISEASE (CVD)	A disease state comprised of disorders ranging from tired, achy legs, varicose veins, skin issues to venous leg ulcers (open and healed).	SIGVARIS GROUP manufactures MEDICAL and CONSUMER HEALTH products to help reduce the risk of, improve symptoms of, and provide treatment for CVD.
DOUBLE-COVERED YARN	Spandex fibers are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort.	Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort.
DEEP VEIN THROMBOSIS (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.	SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening.
GRADUATED COMPRESSION VS. ANTI-EMBOLISM	Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for non-mobile patients.	SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation.
HYPOALLERGENIC	Designed to minimize the chance of an allergic reaction.	SIGVARIS GROUP offers non-latex products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin.
ISO CERTIFIED	SIGVARIS GROUP is ISO (International Organization for Standardization) certified, indicating we have a quality system in place that meets the requirements to manufacture Class II medical devices.	SIGVARIS GROUP follows ISO Standard 13485:2003 for our products manufactured in North and South America.
MICROFIBER	An extremely fine synthetic fiber knitted into textiles with the texture and drape of natural-fiber fabric. This fabric has excellent breathability.	Used to create products with a soft and comfortable feel. Find this fiber in 840 Soft Opaque, 830, 143 & 183 Patterns, 820 Microfiber, and 189 BUSINESS CASUAL.
mmHg (MILLIMETERS OF MERCURY)	Defined as the extra pressure generated by a column of mercury one millimeter high; the same measurement used for blood pressure.	SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement.
MERINO WOOL	Merino is the softest and finest wool available in the world and is inherently itch-free.	SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in 422 Thermoregulating Wool, 421 MERINO OUTDOOR PERFORMANCE SOCKS, 240 Merino Wool, and 152 & 192 ALL-SEASON MERINO WOOL.
SEA ISLAND COTTON	Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness.	SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS.
SENSINNOV	A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day.	This top band is featured in 780 Sheer, 840 Soft Opaque, 710 Patterns, and 820 Microfiber to help keep thigh-high stockings in place all day long.
Supima COTTON <b>SUPIMA</b> WORLD'S FINEST COTTONS	Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fiber.	SIGVARIS GROUP uses this premium cotton yarn in the 230 Cotton, 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS, 146 & 186 CASUAL COTTON, and 230, 142 & 182 Cushioned Cotton products, making them comfortable and ideal for those with sensitive skin.



# SIGVARIS GROUP Dealer Portal

The SIGVARIS GROUP Dealer Portal is your online location for securely ordering SIGVARIS GROUP products and managing your account. Place and track your orders, review order history, contact Customer Care, and much more at [sigvaris-online.com](http://sigvaris-online.com).

## DEALER PORTAL BENEFITS


**24/7** access-order any time!

Three ways to securely place an order online:  
Catalog, Quick Entry or Excel Upload

Access to order history and tracking

Shipping options including drop ship and easy access to an address table for quick shipment of reorders

Check your account balance and print invoices



Email [sigvariscares@sigvaris.com](mailto:sigvariscares@sigvaris.com) and you will receive your active account within one business day!

## SIGVARIS GROUP Education



SIGVARIS GROUP has several new and exciting educational opportunities for 2019 to help you get a leg up on leg health knowledge. These opportunities include online learning, webinars, and certified fitter in-services, as well as an accredited basic or advanced course, in a city near you. Learn more at [sigvariseducation.com](http://sigvariseducation.com).

## Become A SIGVARIS Certified Fitter!

Contact your Territory Manager

or call Customer Care at **1-800-322-7744**, ext. 1

or email [sigvariscares@sigvaris.com](mailto:sigvariscares@sigvaris.com)

SIGVARIS  
GROUP





**ORDERING**

Online: [sigvaris-online.com](http://sigvaris-online.com)

Call: 1-800-322-7744, ext. 1

Fax: 1-800-481-5488

**Customer Service:**

Monday–Friday, 8am–7pm (EST)

**QUALITY STANDARD**

Certified Quality System ISO 13485;  
all products are 100% inspected.

[sigvaris.com](http://sigvaris.com)



**Sigvaris**USA

**SIGVARIS, INC.**  
1119 Highway 74 South  
Peachtree City, GA 30269  
Tel 800-322-7744  
[sigvaris.com](http://sigvaris.com)

**SIGVARIS**  
GROUP

## ORDERING

Online: [sigvaris-online.com](http://sigvaris-online.com)

Call: 1-800-322-7744, ext. 1

Fax: 1-800-481-5488

Customer Care:

Monday–Friday, 8am–7pm (EST)

## QUALITY STANDARD

Peachtree City, GA:

Certified Quality System ISO 13485

FDA Registered

All products are 100% inspected.

SIGVARIS GROUP, PFS<sup>2</sup>, Sensinnov, Compreflex, and Accutab are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2019 SIGVARIS, Inc.

HOLD UP is a trademark of HOLD UP is a trademark of Spec Enterprise LLC.

Lycra is a registered trademark of INVISTA.

Supima is a registered trademark of the Supima Association of America Corporation.

NILIT is a registered trademark of Nilit. All rights reserved.

*drirelease* and FreshGuard are copyright ©2019 *drirelease*. *drirelease* is an Optimer Brands company. All Rights Reserved.

Sim-Slide and Magnide are registered trademarks of Arion Holding B.V.