


lymphology
measuring for custom-made
flat-knit toe caps

effective September 1, 2017

Take the following measurements while patient is weight bearing.

Measuring point A-A1




Length (l) A

Medial length: from base joint of big toe to start of arch.

Length (l) A1

Lateral length: From base joint of little toe to start of arch.

Measuring point A & A1




Circumference (c) A

Measure straight across base joint of little toe and base joint of big toe .

Circumference (c) A1

Measure straight across middle of foot (final length).


Measuring point toes 1-5



Circumference (c) X and Z

Generally measure dimensions X and Z to the millimeter and without tension. The X measurement should not be smaller than the Z measurement, and, where necessary, is adjusted to the Z measurement. Ok to add up to 5mm to cX measurement.

Measuring point Toes 1-5




Length (l XZ)

Shortest length from web X to desired length Z.


Customer name _____ Bill to: _____ Ship to: _____

Account # _____ P.O. # _____

Patient Name _____



Foot		Material	
_____ Left quantity		<input type="checkbox"/> 550	
_____ Right quantity		<input type="checkbox"/> mondi	
Style (Choose One)			
<input type="checkbox"/> Seamless toe cap (as attachment) must be same compression as stocking			
<input type="checkbox"/> Seamed toe cap (as attachment) can be a different compression than stocking.			
<input type="checkbox"/> Individual toe cap			
Standard Colors		Trend Colors*	
<input type="checkbox"/> Caramel	<input type="checkbox"/> Anthracite	<input type="checkbox"/> Magenta	
<input type="checkbox"/> Cashmere	<input type="checkbox"/> Aqua	<input type="checkbox"/> Mint-Green	
<input type="checkbox"/> Black	<input type="checkbox"/> Blue-Jeans	<input type="checkbox"/> Moss-Green	
<input type="checkbox"/> Sand	<input type="checkbox"/> Cherry-Red	<input type="checkbox"/> Navy	
	<input type="checkbox"/> Coral		
mediven 550 Fashion Elements* <i>Must pick 1 color & 1 pattern</i>		mediven 550 Design Elements	
Colors Patterns		<input type="checkbox"/> Stars	
<input type="checkbox"/> Berry	<input type="checkbox"/> Animal	<input type="checkbox"/> Ribs	
<input type="checkbox"/> Brown	<input type="checkbox"/> Crosses	<input type="checkbox"/> Pyramids	
<input type="checkbox"/> Gray	<input type="checkbox"/> Ornaments		
Compression		<input type="checkbox"/> Hallux ease (550 only; seamless toe cap only)	
<input type="checkbox"/> CCL I (18-21 mmHg)			
<input type="checkbox"/> CCL II (23-32 mmHg)			
<input type="checkbox"/> CCL III (34-46 mmHg) (550 only)			
Toe Options		Lymphpad (Optional)	
<input type="checkbox"/> Open Toes		_____ Width (cm)	
<input type="checkbox"/> Closed Toes		_____ Length (cm)	
<input type="checkbox"/> without Small Toe		_____ Location	



*Requires an extra 5 days for production.

medi USA

6481 Franz Warner Parkway
Whitsett, NC 27377

mediusa.com | phone 1.800.633.6334 | fax 1.888.570.4554