Silver and Leg Health

Today's health care providers are facing escalating costs and treatment time for lower extremity patient care. In the United States alone, over 80 million people have venous disease¹. These patients are at increased risk for underlying problems such as leg or foot pain, infections, swelling, skin irritations and ulcers.

The progression of venous disease states usually involves skin breakdown and slow healing times; often leading to infections and/or cellulitis. In some cases these conditions are more severe and chronic.


Many people suffering from venous disease find relief using graduated compression stockings.

Juzo Silver stockings offer the therapeutic benefits of compression, combined with the added benefits of X-Static[®].

X-Static

X-Static[®] fibers are made using 99% pure silver that is bonded to the textile fiber. Today X-Static[®] silver fibers are used in a wide range of products including wound care dressings, sports apparel, military combat uniforms, and footwear.







Juzo USA

P.O. Box 1088

www.juzousa.com

support@juzousa.com

Juzo Canada Ltd. 1100 Burloak Dr. Ste. 300 Burlington, ON L7L 6B2 www.juzo.ca support@juzo.ca

Note: All silver garments are recommended to be removed prior to cardiac defibrillation & MRI scans.


1 American College of Phlebology, 2 American Diabetes Association, 3 Testing performed by Noble Fiber Technologies Advanced Material Development Laboratory. Data on file, 4 Antimicrobial testing performed by North American Science Associates. Data on file, 5 Samson RH, Showalter DP, Stockings and the prevention of recurrent venous ulcers Dermatologic Surgery, 1996 Apr; 22(4):373-6, 6 Department of Health & Human Services, CMS, Program Memorandum AB-03-090, June 20,2003, © X-Static and X-Static - The Silver Fiber are registered trademarks of Noble Fiber Technologies, ® Juzo, Keep Moving, Fibersoft & Signature Collection are registered trademarks of Julius Zorn, Inc., © 2020 - Julius Zorn, Inc





Juzo Silver Garments

For Lymphedema




Safe and Natural

Juzo Silver compression stockings are made using X-Static. Silver is safe, natural, and for centuries has been medically recognized for its antimicrobial properties. Juzo Silver garments do not use chemicals to provide antimicrobial properties.

Ionic Release

The antimicrobial effect of silver is generated by the release of silver ions. A release rate of at least one part per million (PPM) is required to offer effective antimicrobial efficacy. Juzo_® Silver products offer a release rate well above the minimum requirement for antimicrobial activity and solely protect the finished product from microbial growth.




Heat Dissipating

Silver is one of the most thermally conductive elements in the environment. Heat between the skin and garment is dissipated into the environment via the silver fibers. This promotes a more pleasant wearing experience.

Wash Resistant

Juzo Silver is permanently bonded to the textile fibers. Laboratory tests show that Juzo Silver compression stockings maintain high efficacy even after 100 washes^{3, 4}.



*Antimicrobial protects solely the finished garment from microbial growth & odor

Laboratory Tested-Odor free

Juzo Silver compression stocking were tested by North American Science Associates (NAMSA) for antimicrobial activity. The tests show that Juzo Silver compression fabric protected itself from Staphylococcus Aureus in less than 4 hours₄, helping to reduce odor from the garment.

Methicillin-resistant Staphylococcus Aureus (MRSA) is a type of bacteria that is resistant to certain antibiotics. These antibiotics include methicillin and other more common antibiotics such as oxacillin, penicillin and amoxicillin. Laboratory tests show that Juzo Silver products are highly effective at protecting the finished product from MRSA⁴.

Juzo Silver Strong & Expert Silver

Juzo Silver Strong & Expert Silver lymphedema garments feature a dense knit that provides excellent containment. Each garment is custom made to provide an optimal fit.

Dense Knit

The dense knit of Juzo Silver Strong & Silver Expert lymphedema garments provide excellent containment and is ideal for moderate to severe cases of edema and lymphedema.

Made-to-Measure, Seamed for a Precise Fit Juzo Silver Strong & Silver Expert garments are made to an individual's exact measurements. Due to its flat knit design, Juzo Silver Strong & Silver Expert can be made to virtually and shape or size.





Available sizes: • 18-21 mmHq • 23-32 mmHq 34-46 mmHa • 50+ mmHa







Available sizes: 18-21 mmHg 23-32 mmHg 34-46 mmHq

Juzo Silver Dynamic

Juzo Silver Dynamic lymphedema garments provide an optimal combination of therapeutic compression and an inconspicuous appearance.

Seamless

Juzo Silver Dynamic lymphedema garments are circular knitted and provide excellent containment.

Therapeutic Knit

Juzo Silver Dynamic lymphedema garments offer a denser knit than other widely used circular knit garments. This enables Juzo Silver Dynamic garments to provide good containment for mild to moderate cases of edema and lymphedema.

Made-to-Measure for a Precise Fit Juzo Silver Dynamic garments are available made-to-measure for a custom fit.





Available sizes: • 20-30 mmHg • 30-40 mmHg 40-50 mmHa



Juzo Silver Ulcer Liner

The Juzo Silver Ulcer Liner, made with X-Static silver, offers mild compression and can be worn by itself or used as an under layer for a Juzo compression stocking, wrap or bandage.

 Mode 9641 Series Compression Mild Style Knee High Sizing S.L





Safe and naturally antimicrobial compression